

DOCUMENT RESUME

ED 107 968

CE 004 181

AUTHOR McGee, Patricia; And Others
 TITLE An Analysis of the Laboratory Assisting Occupation.
 INSTITUTION Ohio State Dept. of Education, Columbus. Div. of Vocational Education.; Ohio State Univ., Columbus. Trade and Industrial Education Instructional Materials Lab.
 SPONS AGENCY Office of Education (DHEW), Washington, D.C.
 PUB DATE [75]
 NOTE 246p.; For related documents, see CE 004 160-180, CE 004 182-206, CE 004 263-268, and CE 004 425-427
 EDRS PRICE MF-\$0.76 HC-\$12.05 PLUS POSTAGE
 DESCRIPTORS Communication Skills; *Job Analysis; Knowledge Level; *Medical Laboratory Assistants; *Medical Technologists; *Occupational Information; Safety; Skill Analysis; Skill Development; Skilled Occupations; Subprofessionals; *Task Analysis; Task Performance; Work Attitudes

ABSTRACT

The general purpose of the occupational analysis is to provide workable, basic information dealing with the many and varied duties performed in the laboratory assistant occupation. The document opens with a brief introduction followed by a job description. The bulk of the document is presented in table form. Eleven duties are broken down into a number of tasks and for each task a two-page table is presented, showing on the first page: tools, equipment, materials, objects acted upon; performance knowledge (related also to decisions, cues and errors); safety--hazard; and on the second page: science; math--number systems; and communications (performance modes, examples, and skills and concepts). The duties include: performing laboratory techniques, hematology tests, clinical chemistry tests, urinalysis, bacteriology procedures, blood bank and serology procedures, histology procedures, and EKG; operating laboratory equipment; utilizing communication skills; and collecting blood. (BP)

 * Documents acquired by ERIC include many informal unpublished *
 * materials not available from other sources. ERIC makes every effort *
 * to obtain the best copy available. nevertheless, items of marginal *
 * reproducibility are often encountered and this affects the quality *
 * of the microfiche and hardcopy reproductions ERIC makes available *
 * via the ERIC Document Reproduction Service (EDRS). EDRS is not *
 * responsible for the quality of the original document. Reproductions *
 * supplied by EDRS are the best that can be made from the original. *

Occupational Analysis

CE004181

ED107968

LABORATORY ASSISTANT

2

U.S. DEPARTMENT OF HEALTH
EDUCATION & WELFARE
NATIONAL INSTITUTE OF
EDUCATION
THIS DOCUMENT HAS BEEN REPRODUCED EXACTLY AS RECEIVED FROM THE PERSON OR ORGANIZATION ORIGINATING IT. POINTS OF VIEW OR OPINIONS STATED DO NOT NECESSARILY REPRESENT OFFICIAL NATIONAL INSTITUTE OF EDUCATION POSITION OR POLICY.

Instructional Materials Laboratory
Trade and Industrial Education
The Ohio State University

5239

AN ANALYSIS OF THE LABORATORY ASSISTING OCCUPATION

Developed By

Patricia McGee
Instructor, Lab Assisting
Eastland J.V.S.
Columbus, Ohio

Ronald L. Moore
Instructor, Lab Assisting
Eastland J.V.S.
Columbus, Ohio

Gordon E. Pickett
Consultant
Batelle Memorial Institute
Columbus, Ohio

Occupational Analysis
E.P.D.A. Sub Project 73402
June 1, 1973 to December 30, 1974
Director: Tom L. Hindes
Coordinator: William L. Ashley

The Instructional Materials Laboratory
Trade and Industrial Education
The Ohio State University

“The activity which is the subject of this report was supported in whole or in part by the U.S. Office of Education, Department of Health, Education, and Welfare. However, the opinions expressed herein do not reflect the position or policy of the U.S. Office of Education, and no official endorsement by the U.S. Office of Education should be inferred.”

TABLE OF CONTENTS

Foreword	v
Preface	vii
Acknowledgment	ix
Job Description	xi
Duties	
A Performing Laboratory	1
B Operating Laboratory Equipment	29
C Utilizing Communication Skills	69
D Collecting Blood	85
E Performing Hematology Tests	95
F Performing Clinical Chemistry Tests	127
G Performing Urinalysis	165
H Performing Bacteriology Procedures	189
I Performing Blood Bank and Serology Procedures	211
J Performing Histology Procedures	221
K Performing EKG	233

5

FOREWORD

The occupational analysis project was conducted by The Instructional Materials Laboratory, Trade and Industrial Education, The Ohio State University in conjunction with the State Department of Education, Division of Vocational Education pursuant to a grant from the U.S. Office of Education.

The Occupational Analysis project was proposed and conducted to train vocational educators in the techniques of making a comprehensive occupational analysis. Instructors were selected from Agriculture, Business, Distributive, Home Economics and Trade and Industrial Education to gain experience in developing analysis documents for sixty-one different occupations. Representatives from Business, Industry, Medicine, and Education were involved with the vocational instructors in conducting the analysis process.

The project was conducted in three phases. Phase one involved the planning and development of the project strategies. The analysis process was based on sound principles of learning and behavior. Phase two was the identification, selection and orientation of all participants. The training and workshop sessions constituted the third phase. Two-week workshops were held during which teams of vocational instructors conducted an analysis of the occupations in which they had employment experience. The instructors were assisted by both occupational consultants and subject matter specialists.

The project resulted in producing one hundred two trained vocational instructors capable of conducting and assisting in a comprehensive analysis of various occupations. Occupational analysis data were generated for sixty-one occupations. The analysis included a statement of the various tasks performed in each occupation. For each task the following items were identified: tools and equipment; procedural knowledge; safety knowledge; concepts and skills of mathematics, science and communication needed for successful performance in the occupation. The analysis data provided a basis for generating instructional materials, course outlines, student performance objectives, criterion measures as well as identifying specific supporting skills and knowledge in the academic subject areas.

PREFACE

This document describes the duties and tasks performed by a laboratory assistant from a behavioral standpoint, analyzing each as to what the worker does, what mental processes he/she uses and how he/she reacts to the work situation.

The tasks required to accomplish the following are to develop competency in skills and related technical knowledge, to become familiar with chemical and medical terminology, to develop desirable work habits and attitudes necessary to obtain gainful employment, and to work in medical, research and development (r&d), and quality control laboratory occupations.

Although there are no task sheets listing attributes for personal development, this is a very important aspect of a program.

No specific duties were listed except for the medical portion of a program because the tasks are general in nature and provide an opportunity to develop a variety of general skills that are adaptable to a variety of employment situations.

ACKNOWLEDGMENT

We wish to acknowledge the valuable assistance rendered by the following subject matter specialists. They provided input to the vocational instructors in identifying related skills and concepts of each respective subject matter area and served as training assistants in the analysis process during the two-week workshops.

Rollin M. Barber, Psychology
The Ohio State University
Columbus, Ohio

Jodi Beittel, Communications
Columbus, Ohio

Diana L. Buckeye, Mathematics
University of Michigan
Avon Lake, Ohio

Rick Fien, Chemistry
The Ohio State University
Beachwood, Ohio

N.S. Gidwani, Chemistry
Columbus Technical Institute
Columbus, Ohio

Bruce A. Hull, Biology
The Ohio State University
Columbus, Ohio

Donald L. Hyatt, Physics
Worthington High School
Worthington, Ohio

Glenn Mann, Communications
Columbus, Ohio

Jerry McDonald, Physical Sciences
Columbus Technical Institute
Reynoldsburg, Ohio

Colleen Osinski, Psychology
Columbus Technical Institute
Columbus, Ohio

David Porteous, Communications
University of Connecticut
Colchester, Connecticut

James A. Sherlock, Communications
Columbus Technical Institute
Columbus, Ohio

Jim VanArsdall, Mathematics
Worthington High School
Worthington, Ohio

Lillian Yontz, Biology
The Ohio State University
Caldwell, Ohio

Acknowledgment is extended to the following I.M.L. staff members for their role in conducting the workshops; editing, revising, proofing and typing the analyses.

Faith Justice
Sheila Nelson
Marsha Opritza -
Rita Buccilla
Peg Bushelman
Carol Fausnaugh
Mindy Fausnaugh
Rita Hastings
Carol Hicks
Sue Holsinger
Barbara Hughes
Carol Marvin
Patti Nye
Kathy Roediger
Mary Salay

Research Associate
Administrative Assistant
Editorial Consultant
Typist
Typist
Typist
Typist
Typist
Typist
Typist
Typist
Typist
Typist
Typist
Typist

JOB DESCRIPTION

A general laboratory assistant performs laboratory tests, applying practical knowledge of one or more fields to problem solving, according to prescribed standards, to determine chemical and physical characteristics or composition of solid, liquid, or gaseous materials. Using independent judgment and discretion in planning lab work, he/she conducts, tests and makes qualitative and quantitative analyses of materials for purposes such as quality control, process control, product development, determining conformity to specifications and the maintenance of health and safety standards.

A laboratory assistant works under the direction of a biochemist, chemical laboratory chief, analytical chemist, inorganic chemist, organic chemist, physical chemist, metallurgist, pathologist or medical technologist to carry out assigned tasks. He/she learns progressively more comprehensive duties and advances in ability to apply knowledge. He/she assists in supervision, training and development of lower level personnel.

A laboratory assistant prepares chemical solutions according to standard formulas. He/she also sets up, operates and adjusts laboratory equipment, such as ovens, gas cylinders, kilns, vacuum chambers, grinders, agitators, centrifuges, and condensers to prepare material for testing.

In industrial testing the laboratory assistant performs physical tests according to established procedures on dry and liquid substances used as ingredients in adhesives, lubricants, paint, paper and other products for purity, viscosity, density, absorption, melting point and flash point using tension balance, pH meter, and other instruments. He/she tests samples of manufactured products to verify conformance with heat resistance, tensile strength, ductility and other specifications.

A laboratory assistant in the medical laboratory performs routine tests in treatment and diagnosis of disease. He/she also prepares tissue samples for pathologist, takes blood samples and executes such laboratory tests as urinalyses and blood counts, using microscope, spectrophotometer, and similar instruments.

JOB DESCRIPTION CONTINUED

A laboratory assistant records and interprets operating and test data. By personal observation and investigation, he/she determines conformance to established procedures, methods and standards. He/she prepares written test reports, including graphs and charts describing procedures used, results obtained and conclusions reached.

11

Duty A Performing Laboratory Techniques

- 1 Clean glassware
- 2 Separate substances by filtration and decantation
- 3 Separate substances by evaporation and distillation
- 4 Determine density
- 5 Prepare solutions
- 6 Titrate solutions
- 7 Standardize solutions
- 8 Purify by coagulation and sedimentation
- 9 Determine boiling point
- 10 Determine melting point
- 11 Analyze by qualitative methods
- 12 Perform chromatography
- 13 Prepare a dispersion

(TASK STATEMENT) CLEAN GLASSWARE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Dirty glassware Cleaning liquid Soap or detergent Brush Wash basin	Determine nature of contaminant Select appropriate cleaning method	Safety Always know characteristics of cleaning solutions Always wear protective clothing and glasses Evacuate fumes Properly handle glassware Hazard Chemical burn Lacerations Inhalation of noxious fumes Explosion and/or fire
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

(TASK STATEMENT) CLEAN CLASSWARE

SCIENCE	MATH - NUMBER SYSTEMS
<p>Forces acting on a body immersed or floating in a liquid [solubility] Possible chemical reactions of cleaning materials Arrangement of molecules, atoms and ions and the effect on structure and strength of materials Necessary concentrations of solutions Transfer of heat from one body to another</p>	
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Listening</p>	<p><u>EXAMPLES</u></p> <p>Instructions Instructions</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension Logic</p>

(TASK STATEMENT) SEPARATE SUBSTANCES BY FILTRATION AND DECANTATION

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Funnel Filter paper or medium Flask or beaker Suction flask Vacuum source	Filter Porosity of filter medium Decant	Safety Use glassware correctly
<u>DECISIONS</u> Select proper technique	<u>CUES</u> Cloudy filtrates indicate lack of precipitant Heavy solutes can be decanted	<u>ERRORS</u> Repeat the process

(TASK STATEMENT) SEPARATE SUBSTANCES BY FILTRATION AND DECANITATION

SCIENCE		MATH – NUMBER SYSTEMS	
Porosity phenomenology Density Viscosity			
COMMUNICATIONS			
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>	
Viewing	Substance for clarity	Visual analysis, Detail and inference, Describing, Color discrimination	

(TASK STATEMENT) SEPARATE SUBSTANCES BY EVAPORATION AND DISTILLATION

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Evaporation container Heat source Condenser Water source Flasks Fume hood Thermometer Stoppers Adapters Boiling chips Lab record book	Evaporation Simple distillation Fractional distillation	Safety Approach dry point cautiously Use glassware correctly Use boiling chips to avoid hot spots Wear protective devices such as Glasses Hazards Burns Explosions
<u>DECISIONS</u> Determine temperature range in which collection takes place	<u>CUES</u> Observe temperature plateaus and changes therefrom Observe liquid level Be aware of nature of material	<u>ERRORS</u> Repeat process

(TASK STATEMENT) SEPARATE SUBSTANCES BY EVAPORATION AND DISTILLATION

SCIENCE	MATH - NUMBER SYSTEMS
<p>Effect of heating and cooling on state of matter [change of matter from one form to another] Different substances have different boiling points</p>	<p>Measure of temperature [to include Kelvin]</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing Viewing</p>	<p><u>EXAMPLES</u></p> <p>Reference material Lab record book Observe liquid level and temperature changes</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, speed/rate, terminology Penmanship, spelling, terminology Visual analysis, Detail/inference</p>

(TASK STATEMENT) DETERMINE DENSITY

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Lab record book Analytical balance Graduated glassware Bottles or beakers Metric ruler Unknown substance Hydrometer	Measure volume accurately Weigh accurately Calculate density or Displace water and calculate	Safety Use glassware correctly
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

DETERMINE DENSITY

(TASK STATEMENT)

SCIENCE	MATH — NUMBER SYSTEMS
<p>Composition of matter, including protons, neutrons, electrons, atoms, molecules, elements</p> <p>Chemical formulas</p> <p>Periodic table</p>	<p>D = weight/volume (48) in gm/ml (solids, liquids) or gm/l (gases)</p> <p>Basic arithmetic skills</p> <p>Measure of metric weight</p> <p>Metric liquid and dry measures</p> <p>Determination of area and volume of cylinders</p> <p>Determination of area and volume of rectangular, cube and right triangular prisms</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p> <p>Writing</p>	<p><u>EXAMPLES</u></p> <p>Measuring scales</p> <p>Reference materials</p> <p>Record data in lab record book</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension</p> <p>Terminology</p> <p>Penmanship, spelling</p>

(TASK STATEMENT) PREPARE SOLUTIONS

24

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Volumetric glassware Balance, analytical Distilled water Chemicals Reagent bottles Laboratory record book	Measure volumetrically Weigh materials Dilute accurately Label accurately (i.e. - chemical formula, strength of solution, date prepared, name)	Safety Always add solute to solvent Handle glassware properly Know characteristics of chemicals Hazard Fumes Burns Lacerations
<u>DECISIONS</u> Determine nature of solution needed	<u>CUES</u> Desired characteristics of product	<u>ERRORS</u> Undesired properties present

(TASK STATEMENT) PREPARE SOLUTIONS

SCIENCE	MATH - NUMBER SYSTEMS
<p>Composition of matter, including protons, neutrons, electrons, atoms, molecules, and elements</p> <p>Concept of solution</p> <p>Arrangement of molecules, atoms and ions and the effect on structure and strength of materials</p> <p>Periodic table</p> <p>Chemical formulas</p>	<p>Molarity = $\frac{FW(gms)}{1000ml}$</p> <p>Normality = gram equivalent weight/1000ml</p> <p>% = parts/100 (32)(31)</p> <p>Ratio = 59</p> <p>Basic arithmetic skills</p> <p>Measures of metric weight</p> <p>Metric liquid and dry measures</p> <p>Ratio and proportion</p>
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
<p>Reading</p> <p>Writing</p>	<p>Formulas, labels, reference materials</p> <p>Formulas, labels, results in lab record book</p>
SKILLS/CONCEPTS	
<p>Comprehension, detail/inference, terminology</p> <p>Penmanship, spelling, description, usage</p>	

(TASK STATEMENT) TITRATE SOLUTIONS

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED-UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Burette (or pipet) calibrated Standard solutions Indicators Beakers or flasks Unknown samples Laboratory record book	Prepare unknown sample Operate burette Recognize end point Calculate results Use indicators	Safety Use glasswear correctly
<u>DECISIONS</u> If back titration is possible Repetition is necessary	<u>CUES</u> Recognize approach of end point	<u>ERRORS</u> Incorrect determination

(TASK STATEMENT) TITRATE SOLUTIONS

SCIENCE	MATH - NUMBER SYSTEMS
<p>Chemical reaction occurring such as:</p> <ul style="list-style-type: none"> Neutralization Oxidation - reduction Precipitation Complex ion formation Buffer systems Chemical equations 	<p>Addition and subtraction of whole numbers Multiplication and division with whole numbers Multiplication and division of decimal fractions Rounding off decimals and whole numbers Measure of metric weight Metric liquid and dry measures Given an instrument of measure, determine precision, and/or accuracy with respect to relative error, tolerance and significant digits Calculate mean Calculate denotation from mean Ratio and proportion Volume x normality = volume x normality $\% = V \times N \times \text{eq. wt.} / \text{mg} \times 100 (32)$</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Writing Reading Listening</p>	<p><u>EXAMPLES</u></p> <p>Results in record book Directions and/or procedures Directions and/or procedures</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Penmanship, spelling, usage Comprehension, terminology Concentration</p>

STANDARDIZE SOLUTIONS

(TASK STATEMENT)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Calibrated burette Primary standard Indicators Beakers or flasks Solution to be standardized Laboratory record book Analytical balance</p>	<p>Accurately measure primary standard and dilute Use indicators Titrate with burette Recognize end point Calculate normality</p>	<p>Safety Use glassware correctly</p>
<p><u>DECISIONS</u> Back titrate if necessary Repetition is necessary</p>	<p><u>CUES</u> Recognize approach of end point End points should be alike</p>	<p><u>ERRORS</u></p>

(TASK STATEMENT) STANDARDIZE SOLUTIONS

SCIENCE	MATH - NUMBER SYSTEMS
<p>Equivalent weight of titrant = equivalent weight of desired constituent</p> <p>Chemical reactions occurring such as:</p> <p>Neutralization</p> <p>Oxidation - reduction</p> <p>Precipitation</p> <p>Complex ion formation</p> <p>Chemical equations</p> <p>Periodic table</p>	<p>$V \times N(\text{normality}) = \text{mg. primary standard/equivalent weight}$</p> <p>Multiplication and division with whole numbers</p> <p>Multiplication and division of decimal fractions</p> <p>Rounding off decimal and whole numbers</p> <p>Measure of metric weight</p> <p>Liquid and dry measures [metric]</p> <p>Given an instrument of measure, determine precision and/or accuracy with respect to relative error, tolerance, and significant digits</p> <p>Ratio and proportion</p> <p>Calculate mean</p>
PERFORMANCE MODES	COMMUNICATIONS
<p>Writing</p> <p>Reading</p> <p>Listening</p>	<p><u>EXAMPLES</u></p> <p>Results in record book</p> <p>Directions and/or procedures</p> <p>Directions and/or procedures</p> <p>15</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Penmanship, spelling, usage</p> <p>Comprehension, terminology</p> <p>Concentration</p> <p>26</p>

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Beakers or flasks Chemicals - solution, non-solvent	Decant or filter Alter pH	Safety Glassware used correctly Wear protective clothing and glasses Hazards Flammability of materials
<u>DECISIONS</u> Determine ultimate end point	<u>CUES</u> Cloudiness	<u>ERRORS</u> Stopping too soon

TASK STATEMENT) PURIFY BY COAGULATION AND SEDIMENTATION

SCIENCE	MATH -- NUMBER SYSTEMS
Hydrogen ion concentration Solubility data Effect on heating and cooling on state of matter [change of matter from one form to another] Periodic table	Basic arithmetic skills Scientific notation
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
Writing Reading Listening	Results in record book Directions and/or procedures Directions and/or procedures
	17
SKILLS/CONCEPTS	
Penmanship, spelling, usage Comprehension, terminology Concentration	

(TASK STATEMENT) DETERMINE BOILING POINT

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Thermometer Heat source Vacuum source Hydrometer Glassware, calibrated Lab record book Distillation apparatus Reference material Transfer vehicle</p>	<p>Use distillation apparatus Place thermometer correctly and read Apply heat</p>	<p>Safety Use glassware correctly Wear protective devices Thermometers are fragile</p> <p>Hazard Mercuric vapor poisoning Lacerations Burns</p>
<p><u>DECISIONS</u></p> <p>Determine rate of heating</p>	<p><u>CUES</u></p> <p>Formation of bubbles</p>	<p><u>ERRORS</u></p> <p>Damaged equipment or samples</p>

(TASK STATEMENT) DETERMINE BOILING POINT

SCIENCE	MATH - NUMBER SYSTEMS
<p>Vapor pressure Evaporation Atmospheric pressure Gas laws - Boyles, Charles, Gay-Lussac, Gas-law formula Fluids under pressure [incompressibility, transfer of pressure] Ideal gases Effect of heating and cooling on state of matter [change of matter from one form to another] Effect of heating and cooling on expansion of materials [change of dimensions] Argument of molecules, atoms, and ions and the effect on structure and strength of materials Periodic table Chemical formulas</p>	<p>Measures of temperature [including Kelvin] Conversion of atmospheres to mm of Hg, 1 atm = 760 mm Hg</p>
PERFORMANCE MODES	COMMUNICATIONS
<p>Writing Reading Listening</p>	<p><u>EXAMPLES</u> Results in record book Directions and/or procedures Directions and/or procedures</p> <p><u>SKILLS/CONCEPTS</u> Penmanship, spelling, usage Comprehension, terminology Concentration</p> <p>19</p>

(TASK STATEMENT) DETERMINE MELTING POINT

21

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Lab record book Melting point apparatus Thermometer Heat source, controlled Unknown sample Spatula	Read thermometer Recognize transformation from solid to liquid	Safety Handle equipment correctly Wear protective devices Thermometers are fragile Hazard Burns Mercuric vapor poisoning Lacerations
<u>DECISIONS</u> Determine melting point	<u>CUES</u> Recognize initial signs of melting	<u>ERRORS</u> Misread endpoints

(TASK STATEMENT) DETERMINE MELTING POINT

SCIENCE	MATH - NUMBER SYSTEMS
<p>Melting point concepts Effect of heating and cooling on state of matter [change of matter from one form to another] Effect of heating and cooling on expansion of materials [change in dimensions] Arrangement of molecules, atoms, and ions and the effect on structure of strength of materials Periodic table Chemical equations</p>	<p>Measures of temperature (including Kelvin)</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Writing Reading Listening</p>	<p><u>EXAMPLES</u></p> <p>Results in record book Directions and/or procedures Directions and/or procedures</p> <p>21</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Penmanship, spelling, usage Comprehension, terminology Concentration</p>

ANALYZE BY QUALITATIVE METHODS

(TASK STATEMENT)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Chemicals Unknown samples Heat source Glassware Lab record book Fume hood Indicators Procedures sheets Test tubes Reference materials Timer Centrifuge	Chemically separate groups of elements Perform differential chemical reactions Analyze changes occurring	Safety Use glassware correctly Wear protective devices Use fume hood whenever gases are formed Hazard Lacerations Burns Noxious fumes
<u>DECISIONS</u> Select procedure Select equipment and sequence	<u>CUES</u> Color change recognition Formation of precipitates	<u>ERRORS</u> Damaged equipment Inadequate observation

TASK STATEMENT ANALYZE BY QUALITATIVE METHODS

SCIENCE	MATH - NUMBER SYSTEMS
<p>Specific characteristics of elements and element groups</p> <p>Chemical reactions</p> <p>Composition of matter, including protons, neutrons, electrons, atoms, molecules, elements</p> <p>Effect of heating and cooling on state of matter (change of matter from one form to another)</p> <p>Transfer of heat from one body to another</p> <p>Arrangement of molecules, atoms, ions, and the effect on structure and strength of materials</p> <p>Chemical formulas</p> <p>Chemical equations</p> <p>Periodic table</p>	<p>Measures of temperature [including Kelvin]</p> <p>Liquid and dry measures [metric]</p> <p>Read and interpret charts, tables, and/or graphs</p> <p>Measure of time</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p> <p>Writing</p> <p>Viewing</p>	<p><u>EXAMPLES</u></p> <p>Operating manual, reference material, procedures</p> <p>Record data and results in lab record book</p> <p>Color change</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, description of mechanism definition, terminology</p> <p>Penmanship, spelling, classification, description</p> <p>Visual analysis, Color discrimination.</p> <p>23</p>

(TASK STATEMENT) PERFORM CHROMATOGRAPHY

35

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Separation media Material to be separated Solvent Glassware Fume hood Lab record book	Prepare separation media Prepare material to be separated Calibrate media Load columns	Safety Avoid contact with skin Avoid inhalation of fumes Hazard Fainting
<u>DECISIONS</u> Determine type of media to use Determine completion of separation	<u>CUES</u> Frequently check calibration Type of substance to be analyzed	<u>ERRORS</u> Ruin media No separation

(TASK STATEMENT) PERFORM CHROMATOGRAPHY

SCIENCE	MATH - NUMBER SYSTEMS
<p>Nature of solvents Capillary action Osmosis Gravity</p>	<p>Basic arithmetic skills Measure metric length</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Writing Reading Viewing</p>	<p><u>EXAMPLES</u></p> <p>Results in lab record book Separation procedure</p> <p>25</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Penmanship, spelling Comprehension, technical terminology Visual analysis, Detail/inference</p> <p>26</p>

PREPARE A DISPERSION

(TASK STATEMENT)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Dispersion equipment Liquid Material to be dispersed Container Dispersion media Microscope Laboratory record book	Operate dispersion equipment at optimum efficiency Assess dispersion equipment capability Microscopic examination of dispersion	Safety Do not exceed limitations of equipment Use protective devices Hazard Chemical spills
<u>DECISIONS</u> Determine degree of dispersion required Select dispersion equipment Select dispersion media if required	<u>CUES</u> Recognize desired dispersion characteristics Recognize failure of equipment in terms of desired dispersion	<u>ERRORS</u>

(TASK STATEMENT) PREPARE A DISPERSION

SCIENCE	MATH - NUMBER SYSTEMS
<p>Centrifugal forces developed by bodies in rotation (force tending to discharge material from a rotating body) Forces acting on a body immersed or floating in a liquid Resistance of materials to liquid flow</p>	<p>Measure of particular size Measure with the Metric and English system and convert between them Scientific notation</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Writing Reading Viewing</p>	<p><u>EXAMPLES</u></p> <p>Results in lab record book Instructional material Microscopic examination of dispersion characteristics</p> <p>27</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Penmanship, spelling, description Comprehension, description of mechanism logic Visual analysis, Logic, Detail/ inference</p>

Duty B Operating Laboratory Equipment

- 1 Operate bunsen burner
- 2 Measure with calibrated glassware
- 3 Calibrate burette
- 4 Weigh with balances
- 5 Read barometer and thermometer
- 6 Operate pH meter
- 7 Operate centrifuge
- 8 Operate spectrophotometer
- 9 Operate microscope
- 10 Operate ovens
- 11 Operate timing devices
- 12 Operate autoclave
- 13 Operate muffle furnace
- 14 Measure with micrometer
- 15 Install regulators
- 16 Operate incubators and waterbaths
- 17 Operate dispersator
- 18 Operate microtome
- 19 Operate linear slide rule

30

(TASK STATEMENT) OPERATE BUNSEN BURNER

40

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Bunsen burner Fuel source (gas) Tubing Striker	Identify purpose of each part Trace route of gas Adjust gas for proper flame Adjust air intake for proper flame	Safety Keep away from volatile solvents Turn off gas when flame is out Hazard Burns Inhalation of gas fumes
<u>DECISIONS</u> Determine proper adjustment	<u>CUES</u> Flame is hottest directly above inner cone (oxidizing flame) Blue flame is desirable	<u>ERRORS</u> Insufficient heat Fumes, Smoke

(TASK STATEMENT) OPERATE BUNSEN BURNER

SCIENCE		MATH — NUMBER SYSTEMS
<p>Combustion Transfer of energy from one form to another [potential to kinetic]</p>		Ratio and proportion
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Viewing	Size and color of flame	Visual analysis, Color discrimination
	31	

(TASK STATEMENT) MEASURE WITH, CALIBRATED GLASSWARE

42

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Calibrated glassware such as: graded cylinders flasks beakers pipets burettes centrifuge tube tubes Liquid to measure Lab record book	Pour liquids safely Handle glassware Determine meniscus	Safety Use of glassware correctly Wear protective devices Hazard Lacerations
<u>DECISIONS</u> Select proper device	<u>CUES</u> Cleanliness is essential for accurate measurement	<u>ERRORS</u> Inaccuracies Contamination

SCIENCE		MATH - NUMBER SYSTEMS	
Meniscus		<p>Metric system of measurement</p> <p>Given an instrument of measure, to determine precision and/or accuracy with respect to relative error, tolerance and significant digits</p> <p>Liquid and dry measures (metric)</p> <p>Given a coding system, recognize and identify each unit involved by assigning necessary symbols, numerical or literal</p>	
COMMUNICATIONS			
<p><u>PERFORMANCE MODES</u></p> <p>Writing</p>	<p><u>EXAMPLES</u></p> <p>Record results in lab record book</p>	<p><u>SKILLS/CONCEPTS</u></p> <p>Penmanship, spelling, usage</p>	<p>33</p> <p>42</p>

(TASK STATEMENT) CALIBRATE BURETTE

44

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Burette Distilled water Reaker or flask Balance (top loader) Thermometer Graph paper Straight edge Pencil Laboratory record book	Volumetrically measure distilled water Weigh distilled water Record temperature Repeat at 10 ml. intervals Calculate error Draw calibration graph	Safety Use glassware correctly
<u>DECISIONS</u> Determine if repetition is necessary	<u>CUES</u> Clean equipment is essential Interpret results for approximate accuracy and precision	<u>ERRORS</u> Inaccurate results

(TASK STATEMENT) CALIBRATE BURETTE

SCIENCE	MATH - NUMBER SYSTEMS
<p>Relationship between volume and weight Effect of heating and cooling on expansion of materials [change of dimensions]</p>	<p>Basic arithmetic skills Measure of metric weight Measures of temperature [to include Kelvin] Liquid and dry measures [metric] Development of graphs comparing two complimentary sets of figures Read and interpret charts, tables, and/or graphs Calculate mean Deviation from the mean Given an instrument of measure, determine precision, and/or accuracy with respect to relative error tolerance, and significant digits</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Writing Reading Listening</p>	<p><u>EXAMPLES</u></p> <p>Record results in lab record book Prepare graph of results Procedures and/or directions Procedures and/or directions</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Penmanship, spelling, terminology, logic Comprehension, vocabulary Concentration, logic</p> <p>35</p> <p>45</p>

(TASK STATEMENT) WEIGH WITH BALANCES

26

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Balance Weighing paper or container Transfer device (spatula) Laboratory record book Substance to be weighed Set of weights	Clean balance pans Level balance Adjust to zero Determine tare weight Read and interpret scales	Safety Care in handling chemicals Instrument is delicate and requires careful handling
<u>DECISIONS</u> Accuracy required Trip and triple beam balance - weighs to 1 decimal place - accurate to nearest whole number Top loading balance - weighs to 2 decimal places - accuracy to 1 decimal place Analytical balance - weighs to 4 decimal places - accurate to 3 decimal places	<u>CUES</u> Care must be taken when adding or removing substances when pan is released Add weight in sequential order	<u>ERRORS</u> Inaccurate weight

WEIGH WITH BALANCES

(TASK STATEMENT)

SCIENCE	MATH - NUMBER SYSTEMS
Simple levers	Measure of metric weight Liquid and dry measures [metric] Addition and subtraction of whole numbers Addition and subtraction of decimal fractions Given an instrument of measure, determine precision, and/or accuracy with respect to relative error tolerance, and significant digits
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Instruction manual Results in lab record book</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, vocabulary, terminology Spelling, penmanship</p> <p>37</p>

(TASK STATEMENT) READ BAROMETER AND THERMOMETER

48

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Barometer Thermometer Lab record book	Determine whether values should be read in English or Metric system Read and interpret scales on instru- ments	Safety Bendable Hazards Mercury contamination Lacerations
<u>DECISIONS</u> Determine appropriate scale	<u>CUES</u> Column separation	<u>ERRORS</u> Inaccurate reading

(TASK STATEMENT) READ BAROMETER AND THERMOMETER

SCIENCE	MATH -- NUMBER SYSTEMS
<p>Atmospheric pressure Arrangement of molecules, atoms and ions and the effect on structure of strength of materials Gas laws Barometric pressure Absolute and relative humidity</p>	<p>Measures of temperature [to include Kelvin] Measures of length [metric] Liquid and dry measure [metric]</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Writing Reading</p>	<p><u>EXAMPLES</u></p> <p>Results in lab record book Scales on instruments</p> <p>39</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Spelling, penmanship, definition Comprehension, Detail/inference</p>

(TASK STATEMENT) OPERATE pH METER

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
pH meter Buffer solution, standard KCl, saturated Distilled water Electrodes - reference - standard Beakers Wash bottle Unknown sample Lab record book Transfer device Stirring device	Care of electrodes Calibrate meter Correct for temperature Read unknown pH	Safety Electrodes can not be scratched, touched or allowed to dry out Hazards Contamination from chemicals used
<u>DECISIONS</u> Determine proper calibration	<u>CUES</u> Liquid should be in motion during cali- bration pH meter must be warm Read on proper scale	<u>ERRORS</u> Inaccurate reading

SCIENCE	MATH - NUMBER SYSTEMS
<p>Concepts: Hydrogen ion concentration Acid - base theory Ionization potential</p>	<p>Scientific notation Ratio - proportion</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Results on indicator Results in lab record book</p> <p>41</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Description, terminology Penmanship, vocabulary</p>

(TASK STATEMENT) OPERATE CENTRIFUGE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Centrifuge Head Cups Shields Centrifuge tubes Balance Water for balancing Lab record book	Balance tubes on opposite sides of centrifuge Set speed and time Decant liquid	Safety Keep cover closed to avoid flying glass Wear protective glasses Centrifuges must be cleaned, and lubricated regularly to insure operation Hazard Lacerations
<u>DECISIONS</u> Select proper counter balance Select speed	<u>CUES</u> Be sure tubes are balanced Clean cups and shields if breakage occurs	<u>ERRORS</u> Damage or breakage Improper separation

(TASK STATEMENT) OPERATE CENTRIFUGE

SCIENCE	MATH — NUMBER SYSTEMS	
Centripetal forces developed in bodies in rotation (force tending to pull material toward center of rotating body)		
COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Reading Writing	Speed indicator Record in lab record book	Comprehension, terminology, description of mechanism Penmanship, spelling, terminology
	43	

(TASK STATEMENT) OPERATE SPECTROPHOTOMETER

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Spectrophotometer Cuvettes Lab record book Standards	Standardize and calibrate Set wave length Record % T or OD	Safety Spillage must be cleaned to prevent damage to instrument All parts must be sealed and tightened to function properly
<u>DECISIONS</u> Select wave length	<u>CUES</u> Meter must be warm Cuvettes must be free from scratches or fingerprints Bulbs burn out easily Cover must be closed when reading	<u>ERRORS</u> No reading

(TASK STATEMENT) OPERATE SPECTROPHOTOMETER

SCIENCE		MATH -- NUMBER SYSTEMS
<p>Fundamentals of color Beer's law Composition of matter, including protons, neutrons, electrons, atoms, molecules, elements Structure arrangement of molecules, atoms, and ions and the effect on structure and strength of materials Analytical methods</p>		
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
<p>Reading Writing</p>	<p>Instruction manual Results on indicator Results in lab record book</p>	<p>Comprehension, terminology, description of mechanism Penmanship, spelling, logic</p>
	45	55

(TASK STATEMENT) OPERATE MICROSCOPE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Microscope with light source Prepared slide Immersion oil Lens paper Xylo1</p>	<p>Place slide on stage Focus and examine on low power Focus and examine on high power Add oil drop and switch to oil immersion lens</p>	<p>Safety Never focus down on slide With oil, do not touch lens to slide Clean surface with neutral soap and water</p> <p>Hazard Break slide Scratch or break objective</p>
<u>DECISIONS</u> Select objective	<u>CUES</u> Focus using both oculars - use both eyes Do not touch eyepiece with eyelashes Adjust light for best viewing Always use fine adjustment for final focusing Clean eyepieces and objectives before and after use with lens paper Be aware of functions of each part and possible adjustments Higher the magnification, the more light is necessary	<u>ERRORS</u> Poor image Improper focus Damage objective Damage specimen

TASK STATEMENT) OPERATE MICROSCOPE

SCIENCE	MATH -- NUMBER SYSTEMS
<p>Fundamentals of optics Refractive index Types of microscopes</p>	<p>Multiplication and division of whole numbers</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Viewing</p>	<p><u>EXAMPLES</u></p> <p>Microscope parts and usage Image through microscope</p> <p>47</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, description of mechanism terminology Visual analysis, Detail</p>

(TASK STATEMENT) OPERATE OVENS

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Oven Thermometer Material to be dried	Set appropriate temperature Determine length of time necessary for drying	Safety Close oven door Use gloves and/or tongs Hazard Burns
<u>DECISIONS</u> Permit oven to come to equilibrium Set temperature Check thermostat	<u>CUES</u> Function of test Pilot indicator	<u>ERRORS</u> Damage sample

SCIENCE	MATH — NUMBER SYSTEMS
Thermodynamics	
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Operating instructions</p> <p>45</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension</p> <p>55</p>

(TASK STATEMENT) OPERATE TIMING DEVICES

20

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Timer intervals Clock Stop watch	Set timing devices Activate timing devices Read timing devices	Safety Do not drop Do not turn hand backwards Hazard Distraction of timing devices
<u>DECISIONS</u> Select type of timer Check operation	<u>CUES</u> Type of procedure Sound or movement	<u>ERRORS</u> Improper timing Lost sequence

(TASK STATEMENT) OPERATE TIMING DEVICES

SCIENCE	MATH - NUMBER SYSTEMS
	Addition and subtraction of whole numbers Multiplication and division with whole numbers
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Timing device</p> <p>51</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Detail/inference</p>

(TASK STATEMENT) OPERATE AUTOCLAVE

22

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Autoclave Water Material to be autoclaved	Degree of sterilization necessary Evacuate air from autoclave Operate using manual Interpret dials and indicators	Safety Steam burns Hazard Burns
<u>DECISIONS</u> Determine proper water level	<u>CUES</u> Usual conditions are 121° C - 15 to 20 minutes	<u>ERRORS</u> Damaged equipment or material Lost time

OPERATE AUTOCLAVE

TASK STATEMENT)

SCIENCE	MATH -- NUMBER SYSTEMS
<p>Reactions of steam under pressure Conditions necessary for adequate sterilization</p>	<p>Pressure gauge Temperature scale</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Operating instructions manual, dials and indicators</p> <p>53</p>
<p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, terminology, detail/inference</p>	

(TASK STATEMENT) OPERATE MUFFLE FURNACE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Furnace Temperature measuring device Asbestos gloves Long tongs Crucibles Material to be fired Analytical balance</p>	<p>Set appropriate temperature Determine length of time necessary to operate</p>	<p>Safety Use protective devices Use tongs to insert and remove specimens Do not exceed temperature range of furnace Do not handle crucible with hands Hazard Burns</p>
<p><u>DECISIONS</u></p> <p>Permit furnace to equilibrate Determine end point in firing process</p>	<p><u>CUES</u></p> <p>Gages and light cycle</p>	<p><u>ERRORS</u></p> <p>Burn out elements</p>

(TASK STATEMENT) OPERATE NUFFLE FURNACE

SCIENCE	MATH – NUMBER SYSTEMS
Indestructability of energy and matter Effect of heating and cooling on expansion of materials (change of dimensions) Composition of matter, including protons, neutrons, electrons, atoms, molecules, elements Transfer of heat from one body to another Arrangement of molecules, atoms, and ions and the effect on structure and strength of materials Resistance of materials to change in shape [fluidity, elasticity, melting and boiling points]	Measures of temperature [to include Kelvin] Measure of metric weight Liquid and dry measures [metric]
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
Reading	Operating instructions, temperature indicator
	55
	SKILLS/CONCEPTS
	Comprehension, terminology, detail/inference

COMMUNICATIONS

(TASK STATEMENT) MEASURE WITH MICROMETER

66

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Outside micrometer caliper Dial indicator	Hold micrometer Read micrometer scale Proper feel of micrometer for drag	Safety Never measure a rotating shaft Protect form exposure to corrosive Health protective devices Hazard Injury to hand
<u>DECISIONS</u> Degree of accuracy required	<u>CUES</u> Surface condition of micrometer: and material to be measured	<u>ERRORS</u> Damage to micrometer or sample

TASK STATEMENT) MEASURE WITH MICROMETER

SCIENCE	MATH - NUMBER SYSTEMS
Work input, work output, friction and efficiency in simple machines [feel for drag] Arrangement of molecules, atoms and ions and the effect on structure and strength of materials [over tightening]	Measure of length [to include metric] Addition and subtraction of decimal fractions Given an instrument of measure, determine precision and/or accuracy with respect to relative error, tolerance and significant digits Multiplication and division with whole numbers

COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading	Scales and/or dial indicators	Detail/inference, comprehension, technical terminology

57

(TASK STATEMENT) INSTALL REGULATORS

68

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Regulator Bottled gas Wrench Tubing Fittings	Remove cap Secure gas cylinder Purge valve Secure regulator to tank Purge regulator Connect transmission line	Safety Hydrocarbon material Releaf valve Open tank valve all the way Protective devices Hazard Possible explosion
<u>DECISIONS</u> Select proper operating pressure Select proper gauge for different type of gases	<u>CUES</u> Type of regulator Type of gases	<u>ERRORS</u> Damage to equipment

ASK STATEMENT) INSTALL REGULATORS

SCIENCE	MATH — NUMBER SYSTEMS	
Gav - Lussac law Charles law Boyles law	Numbering, system Tank pressure Work pressure	
COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Reading	Instruction manuals	Comprehension
		59

(TASK STATEMENT) OPERATE INCUBATORS AND WATERBATHS

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Incubator or waterbath Material to be incubated Water for waterbath Thermometer	Determine purpose of incubation Determine temperature requirements for particular procedure	Safety Keep water away from electrical equipment Hazard Electrical shocks
<u>DECISIONS</u> Select temperature	<u>CUES</u> Make sure temperature remains constant Keep water level constant Temperature 62°s	<u>ERRORS</u> Damage test media Over heating

WASK STATEMENT) OPERATE INCUBATORS AND WATERBATHS

SCIENCE		MATH - NUMBER SYSTEMS
<p>Optimum growth requirements Optimum color development or reaction completion</p>		
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading	Procedures	Comprehension
	61	

(TASK STATEMENT) OPERATE DISPERSATOR

72

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Dispersing equipment Liquid Material to be dispersed Container Dispersing media Microscope Power source</p>	<p>Set up according to direction Operate according to directions Select and assemble proper equipment</p>	<p>Safety Selection of proper electrical source Hazard Electrical shock</p>
<p><u>DECISIONS</u> Select dispersion equipment Select dispersing media if required Determine purpose and expected results</p>	<p><u>CUES</u> Recognize failure of equipment to perform</p>	<p><u>ERRORS</u> Improper dispersion</p>

SCIENCE	MATH - NUMBER SYSTEMS
<p>Centrifugal forces developed by bodies in rotation 'force tending to discharge material from a rotating body' Forces acting on a body immersed or floating in a liquid Resistance of materials to liquid flow measure of particle size</p>	
PERFORMANCE MODES	COMMUNICATIONS
<p>Reading</p>	<p>Instruction manual</p> <p>63</p> <p>SKILLS/CONCEPTS</p> <p>Comprehension</p>

(TASK STATEMENT) OPERATE MICROTOME

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Microtome Knife Hone Strop Material to be microtomed	Insert sharpened knife Insert mounted specimen Adjust for desired thickness Obtain desired sections Hone and strop knife Adjust knife set screws	Safety Handling blade Hazard Severe lacerations
<u>DECISIONS</u> Determine a satisfactory specimen Select type of cutter	<u>CUES</u> Make adjustment cuts (five) Material to be cut	<u>ERRORS</u> Damage to cutter Damage to specimen Improper cut

TASK STATEMENT) OPERATE MICROTOME

TASK STATEMENT) OPERATE MICROTOME	
SCIENCE	MATH - NUMBER SYSTEMS
Inertia and momentum [bodies at rest and bodies in motion] Resistance of materials to change in shape [bending, twisting, stretching]	Measures of length [metric] Given an instrument of measure, determine precision and/or accuracy with respect to relative error tolerance and significant digits
COMMUNICATIONS	
PERFORMANCE MODES	SKILLS/CONCEPTS
Reading	Comprehension, terminology, detail / inference
65	

(TASK STATEMENT) OPERATE LINEAR SLIDE RULE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Slide rule	Check and adjust scale alignment Determine proper scales Manipulate fine adjustment of the slide and cross hair	Safety Handle with care Hazard Misalignment of scales
<u>DECISIONS</u> Scales required for calculation Select proper index	<u>CUES</u> Standard procedures	<u>ERRORS</u> Improper reading

SCIENCE	MATH - NUMBER SYSTEMS
Theory of scale function	<p>Multiplication and division with whole numbers Multiplication and division of decimal fractions Rounding off decimals and whole numbers Extracting square root Understanding and use of logarithms Scientific notation Use of trigonometric functions in solution of problems involving right triangles</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Instruction manual, scale</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, detail/inference</p> <p>67</p>

Duty C Utilizing Communication Skills

- 1 Read and follow specifications (procedures)
- 2 Prepare table of data
- 3 Prepare graphs
- 4 Make oral presentation
- 5 Record data in laboratory record book
- 6 Write reports
- 7 Utilize reference material

78

(TASK STATEMENT) READ AND FOLLOW SPECIFICATIONS (PROCEDURES)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Specifications Laboratory record book Equipment and material called for in specifications	Read for general information Read for specifics Assemble materials and equipment Follow logical sequential order of specified work Report work in specified terms	Hazard Specified results are determined by following procedures specified
<u>DECISIONS</u> Determine and select proper spec- ifications	<u>CUES</u>	<u>ERRORS</u>

30

(TASK STATEMENT) PREPARE TABLE OF DATA

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Laboratory record book Pencil Columned paper Data to be tabulated</p>	<p>Determine purpose of the table Organize data in logical arrangement Label accurately Enter data</p>	
<p><u>DECISIONS</u></p> <p>Size of tables</p>	<p><u>CUES</u></p> <p>Type of data</p>	<p><u>ERRORS</u></p> <p>Difficult to read</p>

TASK STATEMENT) PREPARE TABLE OF DATA

SCIENCE	MATH - NUMBER SYSTEMS
	<p>Report data is specified mathematical terms Use of numbers (without calculation) Indexing Coding</p>
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
<p>Writing Reading</p>	<p>Prepare table: column heading, footnotes, numerical data Information to be tabulated</p> <p>73</p>
SKILLS/CONCEPTS	
<p>Penmanship, informational reports, format/content, usage Detail/inference, speed/rate</p>	

(TASK STATEMENT) PREPARE GRAPHS

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Graph paper Pencil Straight edge French curve Data to be graphed	Determine purpose for the data and use of graph Apply general terminology pertaining to graphs Prepare graphs using a) linear paper, b) semi-log paper, c) log-log paper	
<u>DECISIONS</u> Select size of graph paper	<u>CUES</u> Type of indicator to be used	<u>ERRORS</u> Poor readability

SCIENCE		MATH - NUMBER SYSTEMS
		Read and interpret charts, tables, and/or graphs Development of graphs comparing two complimentary sets of figures Locate by approximation rational numbers and integers on the number line (sequential ordering)
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing Speaking	Graphs, charts Plot graphs With supervisor	Comprehension, recommendation reports Penmanship, format, usage Punctuation, noise, terminology/ vocabulary
75		83

(TASK STATEMENT) MAKE ORAL PRESENTATION

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Laboratory records	<p>Speak clearly and concisely in "trade", language and corresponding "common", language Analyze and interpret laboratory results orally</p>	
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

(TASK STATEMENT) MAKE ORAL PRESENTATION

SCIENCE		MATH – NUMBER SYSTEMS	
COMMUNICATIONS			
<u>PERFORMANCE MODES</u> Speaking	<u>EXAMPLES</u> Presentation	<u>SKILLS/CONCEPTS</u> Clarity of expression, conciseness, technical vocabulary, organization, diction, logic, work usage, emotional appeal	77
			602

(TASK STATEMENT) RECORD DATA IN LABORATORY RECORD BOOK

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Laboratory record book Pen Data to be recorded	<u>Enter</u> Title Date Purpose Procedure Data Observations Calculations Conclusions	
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

SCIENCE	MATH - NUMBER SYSTEMS
	<p>Basic arithmetic skills Measure of metric length & weight, time & speed, temperature [to include Kelvin], liquid & dry measures [metric] Extracting square root Solution of problems involving numerical & literal algebraic expressions Use of exponents to indicate the power of a number Addition of positive and negative numbers Algebraic subtraction, multiplication and division of numerical and literal terms Manipulation of formula involving three factors Ratio and proportion Determination of area and volume of cylinders Development of graphs comparing two complementav sets of figures Given an instrument of measure, determine precision and/or accuracy with respect to relative error, tolerance, and significant digits; Scientific notation Calculate mean, deviation from mean, and relative average deviation (in % or ppt)</p>
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
Reading Writing	Test results Record in lab book 79
SKILLS/CONCEPTS	
Comprehension, terminology, technical vocabulary, detail/inference Penmanship, precision and conciseness	

(TASK STATEMENT) WRITE REPORTS

TOOLS, EQUIPMENT, MATERIALS,
OBJECTS ACTED UPON

Laboratory record book
Paper
Pen

PERFORMANCE KNOWLEDGE

Organizing factual information to:
Write record of telephone conversation
Write letter report
Write informational (progress) report
Write analytical report
Develop data in a logical manner

SAFETY - HAZARD

Hazard
Protection

DECISIONS

Determine audience level

CUES

Nature of data and use of report

ERRORS

Insufficient information

(TASK STATEMENT) WRITE REPORTS

SCIENCE	MATH - NUMBER SYSTEMS
Inductive reasoning	Organize data
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Lab records Informational and analytical reports, Section of annual reports and other external and internal publications</p>
	<p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension Clarity of expression, proper word usage, grammatical and technical precision, progress report, physical experiment</p>

(TASK STATEMENT) UTILIZE REFERENCE MATERIAL

54

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
References Handbook of chemistry and physics Others	Flexibility in interpreting basis Technical vocabulary	
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

92

Duty D Collecting Blood

- 1 Determine type blood sample needed
- 2 Preform capillary puncture
- 3 Perform venipuncture using syringe or vacutainer
- 4 Prepare blood samples for analysis

73

(TASK STATEMENT) DETERMINE TYPE BLOOD SAMPLE NEEDED

24

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Requisition for blood test	<p>Determine time of drawing Interpret from procedure whether serum, plasma, whole blood are needed Determine amount blood necessary Determine type anticoagulant needed dictated by tests to be performed</p> <p>a) Sodium fluoride (glucose) b) Potassium oxalate (chemistry process) c) Double oxalate (hematology) d) EDTA (hematology) e) Sodium citrate (prothrombin) f) Heparin</p>	
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

(TASK STATEMENT) DETERMINE TYPE BLOOD SAMPLE NEEDED

SCIENCE	MATH -- NUMBER SYSTEMS
<p>Anticoagulant theory Blood clotting mechanism Characteristics of different anticoagulants Composition of blood Treatment necessary to obtain: 1) serum 2) plasma 3) whole blood</p>	
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Requests for laboratory work, procedure</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, technical (medical) terminology</p> <p>87</p>

(TASK STATEMENT) PERFORM CAPILLARY PUNCTURE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Alcohol Lancet Cotton Collection equipment for specific tests	Determine collection site Apply antiseptic Puncture site to obtain free flowing sample. Remove blood sample Stop bleeding	<p>Safety Complete sterility must be maintained Used blades must be discarded safely Technician should not 'stick' self</p> <p>Hazard Contamination of sample Infect patient - unsterile conditions Lacerations from used blades Infection of technician</p>
<p><u>DECISIONS</u></p> <p>Determine where to take sample</p>	<p><u>CUES</u></p> <p>Area used must be free from organisms or edema First drop of blood can not be used Excessive squeezing contaminates sample with tissue fluid Pipets should not touch skin</p>	<p><u>ERRORS</u></p> <p>Contamination Injury</p>

97

(TASK STATEMENT) PERFORM VENIPUNCTURE USING SYRINGE OR VACUTAINER

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Alcohol Cotton or gauze Syringe and needle -or- Shield, vacutainer, needle Tourniquet	Prepare patient Prepare needle and syringe Apply the tourniquet Select vein Apply antiseptic Insert needle Withdraw blood Release tourniquet Withdraw needle Prevent bleeding Transfer blood to proper container	Safety Aseptic conditions must prevail Used equipment must be disposed of properly Tourniquet application must not be prolonged - release before withdrawing needle Needle should not go through vein Artery should not be punctured Assure not to puncture self Hazard Infections of patient Hematoma Prolonged bleeding from lack of pressure applied Lacerations
<u>DECISIONS</u> Select proper location	<u>CUES</u> Arterial blood spurts; venous flows Needle and syringe must be dry Needle must be large enough to prevent trauma Blood must be allowed to flow Patient and technician must be in comfortable position	<u>ERRORS</u> Injury Miss vein or puncture vein

(TASK STATEMENT) PERFORM VENIPUNCTURE USING SYRINGE OR VACUTAINER

TASK STATEMENT PERFORM VERIFICATION USING, SIKING, OR VACUUM		MATH — NUMBER SYSTEMS
SCIENCE	Circulatory system—anatomy and physiology Aseptic conditions Patient concern and treatment	
COMMUNICATIONS		
<u>PERFORMANCE MODES</u> Reading Writing	<u>EXAMPLES</u> Requisitions Labels	<u>SKILLS/CONCEPTS</u> Comprehension, technical (medical) terminology, process report - instructions Accuracy, penmanship
	91	99

(TASK STATEMENT) PREPARE BLOOD SAMPLES FOR ANALYSIS

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Blood sample Anticoagulant Centrifuge Test tubes or vacutainers Pasteur pipet and bulb Refrigerator - freezer	Mix whole blood samples with anti-coagulant Mix samples for plasma, centrifuge and separate plasma from cells Prepare serum samples a) allow blood to clot b) centrifuge c) aspirate serum from clot Label specimens completely Store specimens if possible	Safety Balance centrifuge Use pipets carefully Do not aspirate serum into mouth Do not preserve unless absolutely necessary Hazards Distortion of sample Breakage of sample Lacerations Breakdown of constituents
Decide what analysis to be made	<u>CUES</u> Hemolysis necessitates redrawing of specimen Abnormalities in sample should be noted Clean, dry glassware must be used Most determinations must be performed in fresh specimens. If stored, allow to return to room temperature before use	<u>ERRORS</u> Improper results

(TASK STATEMENT) PREPARE BLOOD SAMPLES FOR ANALYSIS

SCIENCE		MATH — NUMBER SYSTEMS	
Composition of blood Characteristics of blood Appearance of blood - normal and abnormal			
COMMUNICATIONS			
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>	
Reading Writing	Instructions Label specimens 		

Duty E Performing Hematology Tests

- 1 Perform RBC using hemacytometer
- 2 Perform WBC using hemacytometer
- 3 Perform hemoglobin by cyanmethemoglobin method
- 4 Perform microhematocrit
- 5 Prepare blood smear
- 6 Stain blood smear
- 7 Examine blood smear
- 8 Perform erythrocyte sedimentation rate by Wintrobe method
- 9 Perform bleeding time (Duke or Ivy)
- 10 Perform coagulation tests
- 11 Perform prothrombin time
- 12 Perform reticulocyte count
- 13 Perform platelet count
- 14 Calculate indices
- 15 Maintain hematology equipment

(TASK STATEMENT) PERFORM RBC USING HEMACYTOMETER

103

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>RBC diluting pipet (thoma) Aspirator tubing and mouthpiece Blood sample Diluting fluid - Hayem's, Gower's solution Mechanical shaker Hemacytometer Hemacytometer cover glass Alcohol or xylene Microscope Cell counter Pipet cleaning solutions - water, alcohol, acetone Gauze Laboratory record book</p>	<p>Dilute blood sample 1:200 Mix and mount on hemacytometer Count RBC in $1/5 \text{ mm}^2$ on microscope with high dry (40x) Calculate RBC/mm^3</p>	<p>Safety Coverglass can be broken if lens touches Solutions cannot be aspirated into mouth Hazard Scratched or cracked lens Broken coverglass Poisoning from chemicals</p>
<p><u>DECISIONS</u></p> <p>In severe anemia, dilution can be altered</p>	<p><u>CUES</u></p> <p>Repeat until 10% agreement reached Glassware must be clean Improper mounting of specimen causes errors in count</p>	<p><u>ERRORS</u></p>

(TASK STATEMENT) PERFORM RBC USING HEMACYTOMETER

SCIENCE	MATH - NUMBER SYSTEMS
<p>Characteristics and purpose of RBC's Formation of RBC's Decrease = anemia Increase = polycythemia Function of hemacytometer Normal = $4.5 - 5.5 \text{ million/mm}^3$ - men = $4.0 - 5.0 \text{ million/mm}^3$ - women Terminology of RBC and diseases Sources of error Parts of CBC Destruction of RBC Variations from normal</p>	<p>Formula - $\frac{\text{Number of cells counted} \times \text{dilution}}{\text{in } 1/5 \text{ mm}} = \text{RBC/mm}^3$ 0.2 area x 0.1 depth or number of cells x 20,000 Ratio and proportion Basic arithmetic skills</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Directions Record results, labels</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, technical (medical) terminology Penmanship, spelling, accuracy</p> <p>97</p>

(TASK STATEMENT) PERFORM WBC USING HEMACYTOMETER

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
WBC diluting pipet (thoma) Aspirator tubing and mouthpiece Blood sample Diluting fluid - 2% acetic acid or 0.1N HCL Mechanical shaker Hemacytometer Hemacytometer cover glass Alcohol Microscope Cell counter Pipet cleaning solutions Gauze Laboratory record book	Dilute sample 1:20 Mix and mount on hemacytometer Count WBC in four corner square mm with low power (10x) Calculate WBC/mm ³	Safety Aspirate blood and diluting fluid carefully Cover glass can be broken if lens touches Hazard Scratched or cracked lens Broken coverglass Chemical poisoning
<u>DECISIONS</u> Identify area to be used in count	<u>CUES</u> Count cells touching top and left linear Repeat to agreement within 500/mm ³	<u>ERRORS</u> Faulty equipment Faulty technique Cell distribution Nature of sample Small sample size

(TASK STATEMENT) PERFORM WBC USING HEMACYTOMETER

SCIENCE	MATH - NUMBER SYSTEMS
<p>Characteristics and purposes of WBC Formation of WBC's Increase = leukocytosis Decrease = leukopenia Function of reagent and equipment Normal = 5000 - 10,000/mm³ Terminology of WBC and diseases Sources of error Parts of CBC</p>	<p>Number of cells counted = dilution of blood = $\frac{\text{WBC/mm}^3}{\text{Volume (area} \times \text{depth)}}$ or number counted x 50 Ratio and proportion Basic arithmetic skills</p>
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
<p>Reading Writing</p>	<p>Procedures Record results, label</p> <p>99</p>
SKILLS/CONCEPTS	
<p>Comprehension, medical terminology Penmanship, spelling, accuracy</p>	

(TASK STATEMENT) PERFORM HEMOGLOBIN BY CYANMETHAEMOGLOBIN METHOD

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Spectrophotometer Drabkin's reagent (fresh) Cuvettes Pipettes Sahlb pipet Tubing and mouthpiece Blood sample Semi-log graph paper Straight edge Commercial hemoglobin standard Laboratory record book Timer	Standardize photometer using commercial standards. Draw calibration curve Dilute blood 1:251 (5 ml reagent .02 blood) Allow color development Read percent transmittence on photometer Determine concentration from calibration curve	Safety Drabkin's reagent is poisonous Hazard Inhalation of fumes
<u>DECISIONS</u> Determine conditions of equipment and reagents	<u>CUES</u> Accuracy is dependent on technique, equipment, reagent stability, glassware cleanliness Control should be utilized to check technique, equipment, reagents Rinse Sahlb pipet with reagent	<u>ERRORS</u> Improper results

(TASK STATEMENT) PERFORM HEMOGLOBIN BY CYANMETHEMOGLOBIN METHOD

SCIENCE	MATH - NUMBER SYSTEMS
<p>Characteristics and purposes of hemoglobin Anemia = decrease in quality and quantity of RBC and hemoglobin Compounds of hemoglobin Terminology and abbreviation - hemoglobin</p>	<p>Liquid and dry measures [metric] Ratio and proportion Read and interpret charts, tables, and/or graphs Measure with the metric system Given an instrument of measure, determine precision and/or accuracy with respect to relative error, tolerance, and significant digits Measure of time Basic arithmetic skills</p>
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
<p>Reading Writing</p>	<p>Procedure Record results, label, draw calibration graph</p> <p>101</p>
SKILLS/CONCEPTS	
<p>Comprehension, medical terminology Penmanship, spelling, accuracy, format description</p>	

(TASK STATEMENT) PERFORM MICRO-HEMATOCRIT

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Capillary or whole blood sample Heparinized capillary tubes Plasticene or seal-ease (sealing medium) Micro-hematocrit centrifuge Micro-capillary reader Laboratory record book	Fill capillary tube Seal vacant end Centrifuge Read % on reader	Safety Use caution when handling capillary tubes Do not lift cover until rotation has completely stopped Hazard Trauma - lacerations Trauma - to operator
<u>DECISIONS</u> Determine how to take sample	<u>CUES</u> Blood sample must be properly collected and preserved Centrifugation must be standardized	<u>ERRORS</u> Improper results

(TASK STATEMENT) PERFORM MICRO-HEMATOCRIT

SCIENCE		MATH — NUMBER SYSTEMS	
Uses of hemátocrit Comparison of RBC tests RBC diseases and characteristics Normal = 40 - 50% males 37 - 45% females Centrifugal forced developed by bodies in rotation [force tending to discharge material from a rotating body]	Read and interpret charts, tables, and/or graphs Basic arithmetic skills		
COMMUNICATIONS			
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS	
Reading Writing	Procedures Record results, label	Comprehension, medical terminology Penmanship, accuracy	
		103	

(TASK STATEMENT) PREPARE BLOOD SMEAR

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Capillary blood sample Pricleaned glass slides Spreader slide Pencil	Spread drop of blood with spreader slide thinly Air dry smear Label accurately in blood	Hazard Lacerations from slides
DECISIONS Select methods for preparing smear	CUES Blood must not be contaminated with tissue fluid or antiseptics Feather edge indicates good smear - smooth, moderately thin	ERRORS Damage sample Poor quality smear

(TASK STATEMENT) PREPARE BLOOD SMEAR

SCIENCE	MATH — NUMBER SYSTEMS
Purpose of blood smears Terminology of hematology Characteristics of good smears Sources of blood samples	
COMMUNICATIONS	
<u>PERFORMANCE MODES</u> Writing	<u>EXAMPLES</u> Label
<u>SKILLS/CONCEPTS</u> Penmanship, accuracy	
105	

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Well made blood smear Wright's stain Phosphate buffer (pH 6.4) Water Staining rack Timer	Fix smear with Wright's stain Add buffer to stain Wash with water Wipe dye off back of slide Air dry, standing on end	Hazard Stained hands and cloths
DECISIONS Determine standards for accuracy	CUES Accurate timing is essential Metallic sheen should form with buffer pH must be exact Errors occur due to reagents, timing Timing varies with each batch of reagents Rapid drying prevents distortions	ERRORS Incorrect results

(TASK STATEMENT) STAIN BLOOD SMEAR

SCIENCE	MATH — NUMBER SYSTEMS
<p>Wright's stain fixes dead cells Chemical affinity of tissues for stain pH hypotonic - hypertonic</p>	<p>Measure of time</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Procedure</p> <p>107</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, technical (medical) terminology</p>

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Properly stained blood smear</p> <p>Microscope</p> <p>Immersion oil</p> <p>Blood cell calculator</p> <p>Laboratory record book</p>	<p>Evaluate smear quality with low power objective</p> <p>Estimate white count and scan for abnormal cells (low power)</p> <p>Examine RBC morphology (oil immersion)</p> <p>Evaluate platelets (oil)</p> <p>Identify and count 100 WBC's (oil)</p> <p>Repeat abnormalities and % of each type WBC</p>	<p>Hazard</p> <p>Microscope lens can be scratched or broken</p>
<p><u>DECISIONS</u></p> <p>Questionable smears should be reviewed by pathologist</p>	<p><u>CUES</u></p> <p>Number of cells to be counted is dependent on total WBC</p>	<p><u>ERRORS</u></p>

(TASK STATEMENT)

EXAMINT. BLOOD SMEAR

SCIENCE	MATH - NUMBER SYSTEMS
<p>Types of white cells - Neutrophils, Eosinophils, Basophils Lymphocytes, Monocytes Appearance of normal and abnormal or immature RBC's and WBC's Blood cell developmental series Abnormal RBC conditions - anisocytosis, parkilocytosis, hypochromasia, sickle cells, polychromatophilia, basophilic stippling, nucleated red blood cells, target cell Terminology of hematology Objectives of differential Disorders indicated by abnormalities Correlation with other hematology tests</p>	<p>Correct for nucleated RBC Correction = $\frac{\text{uncorrected}}{\text{nm}} \times 3$ Finding a percent of a number and what percent one number is of another Locate by approximation rational numbers and integers on the number line (sequential ordering)</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing Viewing</p>	<p><u>EXAMPLES</u></p> <p>Morphology guides and procedures Record results Blood smear</p> <p>109</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension Penmanship, spelling Visual analysis, Detail/inference, Color discrimination</p>

(TASK STATEMENT) PERFORM ERYTHROCYTE SEDIMENTATION RATE BY WINTROBE METHOD

117

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Whole blood sample Wintrobe tube Pasteur pipet Sedimentation rack Centrifuge Timer Correction chart Laboratory record book	Fill Wintrobe tube accurately Time vertical standing accurately Read tube graduations Centrifuge Correct for anemia	
<u>DECISIONS</u> Select anticoagulant	<u>CUES</u> Anticoagulant used must preserve RBC morphology No hemolysis can be present Tubes must stand perfectly vertical Test must be performed within 2 hours of drawing Accurate timing is essential	<u>ERRORS</u> Improper results Damaged sample

(TASK STATEMENT) PERFORM ERYTHROCYTE SEDIMENTATION RATE BY WINTROBE METHOD

SCIENCE	MATH - NUMBER SYSTEMS
<p>Plasma proteins</p> <p>RBC morphology</p> <p>Increase - infections, and defense mechanisms active</p> <p>Normal - 5 - 20 mm for women</p> <p>5 - 15 mm for men</p> <p>Sources of error</p>	<p>Read sequential scales</p> <p>Measure of time</p>
COMMUNICATIONS	
<u>PERFORMANCE MODES</u>	<u>SKILLS/CONCEPTS</u>
<p>Reading</p> <p>Writing</p>	<p>Comprehension, technical (medical) terminology</p> <p>Penmanship, spelling, accuracy</p>

(TASK STATEMENT) PERFORM BLEEDING TIME (DUKE OR IVY)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
70% alcohol Lancet Stopwatch Blotting paper Cotton Blood pressure cuff Laboratory record book	Perform capillary puncture Time bleeding accurately	Hazard Technician should not stick self
<u>DECISIONS</u> Decide where to puncture	<u>CUES</u> Adequate, standard-sized punctures are essential (3mm depth) Prolonged or shortened results must be repeated Increase shows platelet decrease	<u>ERRORS</u> Improper results

TASK STATEMENT) PERFORM BLEEDING TIME (DUKE OR IVY)

SCIENCE	MATH — NUMBER SYSTEMS	
Clotting mechanism Tissue factors Normal = 1-6 minutes	Measure of time	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedures Results	Comprehension, technical terminology Penmanship, spelling

(TASK STATEMENT) PERFORM COAGULATION TESTS

124

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Test tubes (ID, 8 mm) Blood sample Stopwatch Non-heparinized capillary tubes Lancet Syringe and needle Alcohol and cotton balls Tourniquet Saline Laboratory record book	Describe clot retraction Interpret capillary coagulation time Accurately read Lee-White clotting time	Safety - use glassware correctly Hazard - lacerations from tubes
Determine proper sample appearance	<u>CUES</u> Tissue juice contamination alters results (shortens time) Tests do not differentiate between clotting deficiencies and anticoagulant presence Capillary tests are unreliable Lee-White is influenced by test tube size and surface, temperature, and agitation	<u>ERRORS</u> Poor sample Improper results
<u>DECISIONS</u>		

(TASK STATEMENT) PERFORM COAGULATION TESTS

SCIENCE	MATH — NUMBER SYSTEMS
<p>Blood coagulation theory Normal Lee-White = 15 - 25 minutes Normal clot retraction = Begin 1 hour Complete 18 hours Capillary coagulation time = 2 - 6 minutes Methods of prolonging or retarding coagulation time</p>	<p>Measure of time</p>
PERFORMANCE MODES	COMMUNICATIONS
<p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Procedures Record results</p> <p>115</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology Penmanship, spelling, accuracy</p>

(TASK STATEMENT) PERFORM PROTHROMBIN TIME

123

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Blood sample preserved with sodium oxalate (0.1M) Stock thromboplastin Saline 0.02 M Ca Cl₂ Control solution Centrifuge Test tubes - 8 x 75 mm 37° C incubator Timer Pipets Light source Laboratory record book</p>	<p>Prepare correct 1:10 blood sample Perform and interpret control solutions Observe clot formation accurately</p>	
<p><u>DECISIONS</u> Determine standards of accuracy</p>	<p><u>CUES</u> Must be performed within four hours of drawing specimen Used often for patients on anticoagulant therapy or screening for deficiencies Samples must be performed in duplicate or triplicate and agree within one second Controls must be used to validate reagents and procedure Accurate timing is essential</p>	<p><u>ERRORS</u> Improper timing Incorrect results</p>

(TASK STATEMENT) PERFORM PROTHROMBIN TIME

SCIENCE	MATH — NUMBER SYSTEMS
<p>Blood coagulation theory and mechanism Normal = 11.2 - 14.5 seconds Mechanism of anticoagulant therapy</p>	<p>Dilutions Measure of time Measure of metric volume Basic arithmetic skills</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Record results, label</p> <p>117</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology Penmanship, spelling, accuracy</p> <p>124</p>

(TASK STATEMENT) PERFORM RETICULOCYTE COUNT

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Brilliant cresyl blue or new methylene blue Sodium citrate and sodium chloride or Sodium oxalate and sodium chloride Capillary or venous blood Test tube 8 x 75mm Pipet Timer Slides Microscope Immersion oil Hand counter Aspirator tubing and mouthpiece Laboratory record book</p>	<p>Prepare stain Mix stain and blood accurately Spread mixture on slides Identify number reticulocytes/number of erythrocyte counted under oil immersion</p>	<p>Hazard Lacerations from slides Cracked lens or slide Aspiration of solutions</p>
<p><u>DECISIONS</u> Identify type of blood cell</p>	<p><u>CUES</u> 2000 RBC's are counted - 1000 on each of 2 slides Count in medium thin portion of slides Precipitated stain must not be confused with reticulocytes Stain should be filtered immediately before use Allowable difference between slides = 5-7 cells</p>	<p><u>ERRORS</u> Count wrong kind of cell Inaccurate count</p>

(TASK STATEMENT) PERFORM RETICULOCYTE COUNT

SCIENCE	MATH - NUMBER SYSTEMS
<p>RBC maturation series RBC morphology and abnormalities Normal = 0.8 - 1.5% retics Formation of RBC Physiology of anemias Supravital staining Uses of reticulocyte count</p>	<p><u>Number reticulocytes counted</u> x 100 = % retics Number erythrocyte counted Measure of metric volume Basic arithmetic skills</p>
<p style="text-align: center;">COMMUNICATIONS</p>	
<p><u>PERFORMANCE MODES</u></p> <p>Writing Reading Viewing</p>	<p><u>EXAMPLES</u></p> <p>Record results, label's Slide</p> <p style="text-align: right;">119</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Penmanship, spelling, accuracy Comprehension, medical terminology Visual analysis, Color discrimination</p>

(TASK STATEMENT) PERFORM PLATELET COUNT

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Rees-Ecker diluting fluid RBC diluting pipet Capillary or venous blood Pipet shaker Blood smears Hemocytometer Gauze Petri plate contained moistened gauze Microscope Hand counter Aspirator tubing and mouthpiece Laboratory record book	Dilute blood with Rees-Ecker Mount specimen or hemocytometer Allow adequate settling Identify and count platelets using high dry in center square millimeter Calculate results	Hazard Aspiration of solution Broken cover glass Scratched or broken lens
<u>DECISIONS</u> Identify type of cells	<u>CUES</u> Platelets are small, attach easily to glassware, clump easily Diluent must be stored in refrigerator and filtered prior to use Glassware must be very clean Counts must be done in duplicate and must agree within 20,000 cells	<u>ERRORS</u> Count wrong type of cell Inaccurate count

(TASK STATEMENT) PERFORM PLATELET COUNT

SCIENCE	MATH — NUMBER SYSTEMS
<p>Blood coagulation theory Formation of platelets Normal = 170,000 - 400,000 cells/mm³ Terminology of abnormal platelet quantity Use of platelet count</p>	<p>Number of platelets counted x 1mm³ x 200 (dilution of blood) 0.1 mm³ Platelets/mm³ or number counted x 200 Measure of metric volume Basic arithmetic skills</p>
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
<p>Reading Writing Viewing</p>	<p><u>SKILLS/CONCEPTS</u> Comprehension, medical terminology Penmanship, spelling, accuracy Visual analysis, Color discrimination</p> <p>Record results, label Slide</p> <p>121</p>

(TASK STATEMENT) CALCULATE INDICES

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Results of RBC Hematocrit Hemoglobin Anemia classifier (blood constants calculator) Laboratory record book	Calculate MCV Calculate MCH Calculate MCHC	129
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

(TASK STATEMENT) CALCULATE INDICES

SCIENCE	MATH - NUMBER SYSTEMS
<p>Comparison of RBC, Hematocrite, Hemoglobin Normal blood value range Values of indices Theory of indices</p>	<p> $MCV = \frac{\text{Hematocrit} \times 10}{\text{RBC (millions)}} \text{ (cu. microns)}$ $MCH = \frac{\text{Hemoglobin} \times 10}{\text{RBC (millions)}} \text{ (uug)}$ $MCHC = \frac{\text{Hemoglobin} \times 100}{\text{Hematocrit} (\%)} \text{ Basic arithmetic skills}$ </p>

COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
<p>Reading Writing</p>	<p>Procedures, classifier Record and report results</p> <p>123</p>	<p>Comprehension, medical terminology, description of mechanism Penmanship, spelling, progress report</p>

(TASK STATEMENT)

MAINTAIN HEMATOLOGY EQUIPMENT

131

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Water Alcohol Acetone Xylol Gauze Vacuum Pipet washer Cleaning jars Aspirator Blood pipet cleaning wires	Clean diluting pipets Clean hemocytometer Clean slides	Hazards Lacerations Burns
<u>DECISIONS</u> Determine proper conditions	<u>CUES</u> Contamination Failure to function	<u>ERRORS</u> Inoperative equipment

132

Duty F Performing Clinical Chemistry Tests

- 1 Perform glucose test by Nelson-Samogyi method
- 2 Perform glucose by Folin-Wu method
- 3 Perform urea nitrogen test by nesslerization
- 4 Perform uric acid test
- 5 Perform creatinine test (Folin method)
- 6 Perform total protein, albumin, globulin test (Biuret method) (TP/AG)
- 7 Perform cholesterol test
- 8 Perform chloride test (Scholes and Scholes)
- 9 Perform carbon dioxide (CO₂) determination (Van Slyke)
- 10 Perform sodium and potassium determination (Flame photometer)
- 11 Perform calcium test (Clark-Callip)
- 12 Perform inorganic phosphorous test
- 13 Perform amylase test (Samogyi method)
- 14 Perform SGOT, SGPT, LNH tests (Sigma)
- 15 Perform alkaline and acid phosphatase tests
- 16 Perform VDB (Van Den Berge) test (Mallay and Foelyn)
- 17 Perform icteues index, thymol turbidity, cephalin cholesterol flocculation
- 18 Perform BSP (Bromsulphonphthalien) test

(TASK STATEMENT) PERFORM GLUCOSE TEST BY NELSON - SAMOGYI METHOD

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Sample of serum, plasma, whole blood urine or CSF Test tubes Pipet Distilled water Barium hydroxide Zinc sulfate Centrifuge Folin-Wu sugar tubes Alkaline copper tartrate Boiling water bath Cold water bath Arsenomol update color reagent Glucose standard Spectrophotometer Cuvettes Graph paper Timer Heat source Controls Lab record book</p>	<p>Prepare protein - free filtrate Place samples in Folin-Wu tubes Perform color reactions and dilute Calibration spectrophotometer Read % T values Determine concentration from graph with mg %</p>	<p>Safety Reliability and accuracy of results Proper glassware handling Chemical handling Sources of error Hazards Patient's well-being hangs in balance Burns - heat plus chemical Lacerations</p>
<p><u>DECISIONS</u></p> <p>Set proper wavelength Check quality control</p>	<p><u>CUES</u></p> <p>Accurate timing is essential Reagents must be fresh Inadequate mixing</p>	<p><u>ERRORS</u></p> <p>Incorrect results</p>

(TASK STATEMENT) PERFORM GLUCOSE TEST BY NELSON-SAMOGYI METHOD

SCIENCE	MATH -- NUMBER SYSTEMS
Digestion and metabolism of carbohydrates Chemistry of carbohydrates Blood sugar homeostasis Kidney, liver, pancreas function Glucose tolerance theory Normal values and glucose curves Protein free filtrates Proper specimen collection times plus preservation Chemical reactions of test procedures Diabetes mellitus Renal threshold	Measure of metric volume and weight Read graph Dilutions Ratio and proportions Liquid and dry measures (metrics) Measure of time Basic arithmetic skills
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
Reading Writing	Procedure Report results, labels, draw graphs
SKILLS/CONCEPTS Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage	

(TASK STATEMENT) PERFORM GLUCOSE BY FOLIN-WU METHOD

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Test tubes Flashes Whole blood Sulfuric acid Sodium tungstate Funnel Filter paper Folin-Wu blood sugar tubes Distilled water Alkaline copper tartrate Boiling water bath Heat source Cold water bath Molybdic acid Spectrophotometer Pipets Cuvettes Graph paper Timer Lab record book Standard glucose	Prepare Folin-Wu filtrate Place samples in Folin-Wu tubes Perform color reaction and dilute Calibrate spectrophotometer Read % T values Determine mg% from graph	Safety Reliability and accuracy of results Proper glassware handling Chemical handling Sources of error Hazards Patients well being Burns - heat and chemical Lacerations
Select procedure	Set proper wavelength Accurate timing is essential Check quality control Inadequate mixing	<u>ERRORS</u> Improper results

137

(TASK STATEMENT) PERFORM UREA NITROGEN TEST BY NESSLERIZATION

138

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Whole blood sample Pipet Test tubes Water bath or incubator Timer Corks Filter paper Funnels Cuvettes Spectrophotometer Nitrogen standard [$(\text{NH}_4)_2\text{SO}_4$] Controls Graph paper Urease solution Urea solution Sulfuric acid Sodium tungstate Nessler's reagent Distilled water Lab record book	Incubate blood with urease Prepare protein free filtrates Develop color Calibrate spectrophotometer Read % T values Determine concentration from graph in mg %	Safety -- Reliability and accuracy of results Glassware handling Chemical handling Sources of error Hazard Patient's well-being Lacerations Chemical burns
Select time and temperature	<u>CUES</u> Time and temperature must be accurate Wavelength setting Check quality control Make fresh reagent	<u>ERRORS</u> Improper results

(TASK STATEMENT) PERFORM URFA NITROGEN TEST BY NESSLERIZATION

SCIENCE	MATH - NUMBER SYSTEMS
<p>Kidney function Digestion and metabolism of proteins Kidney diseases Nitrogen substances in body Liver function Chemical reactions of test procedure Collection and preservation of specimen Enzyme reactions Normal values</p>	<p>Measure of metrics volume and weight Read graph Dilutions Measure of temperature Ratio and proportions Measure of time Basic arithmetic skills</p>
COMMUNICATIONS	
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>
<p>Reading Writing</p>	<p>Procedure Report results, labels, draw graphs</p>
<p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology, description of mechanism Spelling, penmanship, accuracy, progress report, usage</p>	
<p>133</p>	

(TASK STATEMENT) PERFORM URIC ACID TEST

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Serum or plasma Test tube Pipets Sulfuric acid Sodium tungstate Centrifuge Distilled water Sodium carbonate Phosphotungstic acid reagent Spectrophotometer Cuvettes Timer Uric acid standard Graph paper Lab record book Controls	Prepare Folin-Wu filtrate Develop color Calibrate spectrophotometer Read % T Determine concentration from graph in mg %	Safety Reliability and accuracy of results Classware handling Chemical handling Hazard Patient's well-being Lacerations Chemical burns
<u>DECISIONS</u> Select proper wave lengths	<u>CUES</u> Procedure specifications	<u>ERRORS</u> Inaccurate results

TASK STATEMENT) PERFORM URIC ACID TEST

SCIENCE		MATH — NUMBER SYSTEMS
Kidney function Digestion and metabolism of proteins Nitrogen substances in body Chemical reactions of test procedure Normal values Kidney diseases Liver function	Measure of metric volume Ratios and proportions Read graph Dilutions Measure of time Basic arithmetic skills	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedure Report results, labels, draw graphs	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage

(TASK STATEMENT) PERFORM CREATININE TEST (FOLIN METHOD)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Serum or plasma Water Sodium tungstate Sulfuric acid Test tubes Pipets Centrifuge Cuvette Spectrophotometer Timer Creatinine standard Graph paper Lab record book Controls Picric acid Sodium hydroxide	Prepare Folin-Wu filtrate Develop color Calibrate spectrophotometer Read % T values Determine concentration from graph in mg %	Safety Reliability and accuracy of results Glassware handling Chemical handling Picric acid is volatile when hot Hazard Patient's well-being Lacerations Chemical burns Picric acid stains
<u>DECISIONS</u> Determine if test is complete	<u>CUES</u> Check quality control Check wavelength Fresh reagents Procedure specification	<u>ERRORS</u> Improper results

(TASK STATEMENT) PERFORM CREATININE TEST (FOLIN METHOD)

SCIENCE		MATH — NUMBER SYSTEMS
Chemical reactions of test Jaffe reaction Kidney function Digestion and metabolism of proteins Nitrogen substances in body Nitrogen substances in body Kidney diseases Normal values		Ratios and proportions Measure of metric volume Read graph Dilutions Measure of time Basic arithmetic skills
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedure Report results, labels, draw graphs	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage
	137	143

144

(TASK STATEMENT) PERFORM TOTAL PROTEIN, ALBUMIN, GLOBULIN TEST (BIURET METHOD) (TP/AG)

SCIENCE	MATH — NUMBER SYSTEMS	
Digestion and metabolism of protein Protein components and functions Nitrogen substances in body Normal values Chemical reactions of test Significance of findings	Globulin = total protein - albumin A/C = albumin/globulin Measure of metric volume Ratio and proportions Read graph Dilutions Measure of time Basic arithmetic skills	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedure Report results, labels, draw graphs	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage

(TASK STATEMENT) PERFORM CHOLESTEROL TEST

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Serum Volumetrics flask Funnel Filter paper Stoppers Acetone - alcohol Deaker Pipets Heat source Glacial acetic acid Ferric chloride Cuvettes Spectrophotometer Cholesterol standard Timer Lab record book Controls Graph paper	Prepare filtrate Evaporate to dryness Develop color Calibrate spectrophotometer Read % T values Determine concentration from graph in mg %	Reliability and accuracy of results Glassware handling Chemical handling Alcohol-acetone is flammable Hazard Patient's well-being Lacerations Chemical burns Fire
Determine endpoint	<u>CUES</u> Dry equipment Do not scorch residue Evaporation occurs easily Check quality control Check wavelength	<u>ERRORS</u> Improper results

(TASK STATEMENT) PERFORM CHOLESTEROL TEST

SCIENCE	MATH — NUMBER SYSTEMS	
Digestion and metabolism of fats Interpretation of results Chemical reaction of tests Kidney function and diseases Liver function and diseases Types of lipids and functions	Measure of metric volume Ratios and proportions Read graph Dilutions Measure of time Basic arithmetic skills	
COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Reading Writing	Procedures Report results, labels, draw graphs	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage
	141	

(TASK STATEMENT) PERFORM CHLORIDE TEST (SCHOLES AND SCHOLFS)

148

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Serum Water Test tube Sulfuric acid Sodium tungstate Funnel Filter paper Flask, erlenmeyer Diphenyl carbonzone Mercuric nitrate Standard sodium chloride Burette Lab record book Controls</p>	<p>Prepare Folin-Wu filtrate Titrate with mercuric nitrate Calculate results in mg% and/or mg/l Titrate standard NaCl</p>	<p>Safety Reliability and accuracy of results Glassware handling Chemical handling Sources of error Hazard Patient's well-being Lacerations Chemical burns</p>
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> Check quality control</p>	<p><u>ERRORS</u> Incorrect results</p>

148

(TASK STATEMENT) PERFORM CHLORIDE TEST (SCHOLFS AND SCHOLFS)

SCIENCE	MATH - NUMBER SYSTEMS
<p>Electrolyte balance in body Water balance Chemical reactions of tests Titration procedures Normal values Endocrine system - adrenals</p>	<p>Conversion of milliequivalents - mg% $Ml\ Hg(NO_3)_2 \times 100/A = mg/liter$ $Ml\ Hg(NO_3)_2 \times 100/A \times 5.85 = mg/100ml$ $A = ml$ to titrate standard NaCl Measure of metric volume Ratio and proportion Dilution Basic arithmetic skills</p>
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
<p>Reading Writing</p>	<p>Procedures Report results, labels, draw graphs</p> <p>143</p> <p>SKILLS/CONCEPTS Comprehension, medical terminology, description of mechanism, Penmanship, spelling, accuracy, progress report, usage</p>

(TASK STATEMENT) PERFORM CARBON DIOXIDE (CO₂) DETERMINATION (VAN SLIKE)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Serum Van Slyke apparatus Mercury Distilled water Caprylic alcohol Lactic or sulfuric acid Pipet Lab record book</p>	<p>Prepare apparatus Add reagents and serum Release CO₂ Measure volume Calculate volume % at STP conditions</p>	<p>Safety Mercury Glassware handling Hazard Lacerations Chemical burns</p>
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> Measure temperature and pressure Clear of water and air</p>	<p><u>ERRORS</u> Incorrect results</p>

(TASK STATEMENT) PERFORM CARBON DIOXIDE (CO₂) DETERMINATION (VAN SLYKE)

SCIENCE	MATH — NUMBER SYSTEMS	
Electrolyte balance in body Acidosis and alkalosis Gas laws Chemical reactions of test Manometric determinations Normal values Water balance Endocrine system - adrenals	Metric measure of volume Connect to STP conditions Measure of temperature Measure barometric pressure Basic arithmetic skills	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedures Report results, labels, draw graphs	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage
	145	153

(TASK STATEMENT) PERFORM SODIUM AND POTASSIUM DETERMINATION (FLAME PHOTOMETER)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Serum Distilled water Buffer solutions Standard solutions Volumetric flask Pipets Gas source - propane Air source Controls</p>	<p>Set-up photometer Calibrate Read Na and K values</p>	<p>Safety Reliability and accuracy of results Flow meter Open flame Hazard Lacerations Burns</p>
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> Check quality control Contamination of buffers and standards</p>	<p><u>ERRORS</u> Incorrect results</p>

(TASK STATEMENT) PERFORM SODIUM AND POTASSIUM DETERMINATION (FLAME PHOTOMETER)

SCIENCE	MATH - NUMBER SYSTEMS	
Electrolyte balance in body Buffer systems Water balance Endocrine system - adrenals	Mg. % - meg/l Measure of metric volume	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedures Report results, labels, draw graphs	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage
	147	153

(TASK STATEMENT) PERFORM CALCIUM TEST (CLARK-CALLIP)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Serum Distilled water Ammonium oxalate Graduated centrifuge tube Centrifuge Filter paper Ammonium hydroxide Sulfuric acid Pipets Boiling water bath Heat source Potassium permanganate Timer Lab record book Controls</p>	<p>Precipitate oxalate Centrifuge and drain Heat Titrate with KMnO_4 Titrate blank Calculate mg % Ca</p>	<p>Safety Sources of error Reliability and accuracy of results Glassware handling Chemical handling Hazard Patient's well-being Lacerations Stains Burns</p>
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> Reciprocal relation between calcium and phosphorus Check quality control</p>	<p><u>ERRORS</u> Incorrect results</p>

(TASK STATEMENT) PERFORM CALCIUM TEST (CLARK - CALLIP)

SCIENCE	MATH - NUMBER SYSTEMS	
Electrolyte balance Chemical reactions of test Normal values Titration mechanism Endocrine system - parathyroid, thyroid Vitamin metabolism	1 ml $\text{KMNO}_4=0.2$ mg Ca ml KMNO_4 needed - ml KMNO_4 for blank x 0.2 x 100/2 = mg Ca/100ml Metric measure of volume Basic arithmetic skills Measure of temperature Measure of time	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedures Report results, labels, draw graphs	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage

(TASK STATEMENT) PERFORM INORGANIC PHOSPHOROUS TEST

856

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Serum Graduated centrifuge tube Pipets Trichloroacetic acid Funnel Filter paper Test tubes Molybdate II Amino naphthol sulfuric acid (ANS) Distilled water Standard KH_2PO_4 Spectrophotometer Cuvettes Controls Lab record book Timer Graph paper	Prepare filtrate Develop color Calibrate spectrophotometer Read % T value Determine mg % from graph	Safety Reliability and accuracy of results Proper glassware handling Chemical handling Sources of error Hazard Patient's well-being Burns Lacerations
Determine standards for accuracy	Check quality control Set proper wavelength Avoid hemolysis	Incorrect results
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

(TASK STATEMENT) PERFORM INORGANIC PHOSPHORUS TEST

SCIENCE	MATH - NUMBER SYSTEMS
<p>Electrolyte balance Chemical reaction of test Normal values Endocrine system - parathyroid, thyroid Vitamin metabolism Vitamin deficiency diseases</p>	<p>Measure of metric volume Dilutions Basic arithmetic skills Ratios and proportions Liquid and dry measures Measure of time Read graph</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Procedures Report results, labels, draw graphs</p> <p>151</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage</p> <p>157</p>

(TASK STATEMENT) PERFORM AMYLASE TEST (SAMOCYI METHOD)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Starch paste Serum Test tubes Waterbath Pipets Sulfuric acid Sodium tungstate Reagents for glucose test (Folin-Wu) Lab record book Glucose calibration graph Timer</p>	<p>Incubate serum and substrate Prepare filtrate Perform glucose test Read % T on spectrophotometer Read value from glucose calibration graph</p>	<p>Safety Reliability and accuracy of results Glassware handling Chemical handling Sources of error Hazard Patient's well-being Lacerations Chemical burns</p>
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> Must be performed soon after specimen obtained Temperature, pH, timing critical</p>	<p><u>ERRORS</u> Incorrect results</p>

(TASK STATEMENT) PERFORM AMYLASE TEST (SAMOGYI METHOD)

SCIENCE	MATH - NUMBER SYSTEMS	
Carbohydrate digestion and metabolism Function of pancreas Chemical reactions of test Enzyme function and activity Pancreas pathology and diseases Liver function Amylase function	1 amylase unit = 1 mg % sugar Measure of metric volume and weight Read graphs Dilutions Measure of temperature Ratios and proportions Measure of time Basic arithmetic skills	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedure Report results, label, draw graphs	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage
153		159

(TASK STATEMENT) PERFORM SGOT, SGPT, LDH TESTS (SIGMA)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Serum Substrate (Sigma) Waterbath Centrifuge tubes, graduated Timer Color reagent (Sigma) Sodium hydroxide Spectrophotometer Control serum Standard solution Graph paper Cuvettes Lab record book</p>	<p>Incubate serum and substrate Develop color Calibrate spectrophotometer Read % T value Determine concentration from graph in Sigma-Frankel units</p>	<p>Safety Reliability and accuracy of results Glassware handling Chemical handling Sources of error Hazard Patient's well-being Lacerations Chemical burns</p>
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> Time and temperature critical Check wavelength Check quality control</p>	<p><u>ERRORS</u> Incorrect results</p>

(TASK STATEMENT) PERFORM SBOT, SCPT, LDH TESTS (SIGMA)

TASK STATEMENT)		PERFORM SBOT, SGPT, LDH TESTS (SIGMA)	
SCIENCE		MATH -- NUMBER SYSTEMS	
Chemical reactions of tests Enzyme function and activity Heart function Myocardial infarctions Liver function Pulmonary infarctions Preservation of samples Enzyme curves	Measure of metric volume and weight Read graphs Dilutions Measure of temperature Ratios and proportions Measure of time Basic arithmetic skills		
COMMUNICATIONS			
PERFORMANCE MODES		EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedure Report results, label, draw graphs	155	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage

(TASK STATEMENT) PERFORM ALKALINE AND ACID PHOSPHATASE TESTS

0.02

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Serum Alkaline or acid substrate Water bath or incubator Distilled water Trichloroacetic acid Funnel Filter paper Graduated centrifuge tube Molybdate II solution ANS Pipets Spectrophotometer Cuvettes Graph paper Standard KH_2PO_4 Controls Lab record book</p>	<p>Prepare filtrate Incubate Develop color Calibrate spectrophotometer Read % T value Determine concentration from graph in B units Perform inorganic phosphorus test</p>	<p>Safety Reliability and accuracy of results Glassware handling Chemical handling Sources of error Hazard Patient's well-being Lacerations Chemical burns</p>
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> Time, substrate, temperature, pH critical Check quality control Check wavelength</p>	<p><u>ERRORS</u> Incorrect results</p>

(TASK STATEMENT) PERFORM ALKALINE AND ACID PHOSPHATASE TESTS

SCIENCE	MATH - NUMBER SYSTEMS
<p>Function of enzymes and activity Normal values Chemical reactions of tests Optimum conditions for enzyme tests Prostatic carcinoma Bone diseases</p>	<p>Phosphatase - inorganic phosphorus = Alkaline or acid phosphatase Measure of temperature 1 B unit = 1 mg phosphorus/100 ml serum/ 1 hour Measure of metric volume and weight Read graph Dilutions Measure of temperature Ratios and proportions Measure of time Basic arithmetic skills</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Procedure Report results, label, draw graphs</p> <p>157</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage</p> <p>163</p>

(TASK STATEMENT) PERFORM VDB (VAN DEN BERGLE) TEST (MALLAY AND EOELYN)

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Serum Large test tubes Distilled water Methyl alcohol Hydrochloric acid Sulfanilic acid Sodium nitrite Timer Spectrophotometer Cuvettes Graph paper Standard bilirubin Lab record book</p>	<p>Set up color reaction Calibrate spectrophotometer Read direct reaction Read indirect reaction Determine mg% from graph</p>	<p>Safety Reliability and accuracy of results Proper glassware handling Chemical handling Sources of error Hazard Patient's well-being Burns Lacerations</p>
<p><u>DECISIONS</u></p> <p>Determine standards for accuracy</p>	<p><u>CUES</u></p> <p>Reagents must be fresh Timing is critical Avoid hemolysis</p>	<p><u>ERRORS</u></p> <p>Incorrect results</p>

(TASK STATEMENT) PERFORM VDB (VAN DEN BERGLE) TEST (MALLAY AND EGELYN)

SCIENCE	MATH - NUMBER SYSTEMS	
Digestion and metabolism of food Liver function Red cell breakdown Liver diseases Chemical reactions of test Normal values Bile formation and function	Measure of metric volume Dilutions Basic arithmetic skills Ratios and proportions Measure of metric weight Measure of time Read graph	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedure Report results, label, draw graph	Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage

(TASK STATEMENT) PERFORM ICTEUES INDEX, THYMOL TURBIDITY, CEPHALIN CHOLESTEROL FLOCCULATION

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Serum A Test tubes Sodium citrate Standard potassium dichromate Spectrophotometer Cuvettes Pipettes B Test tubes Buffered thymol Pipettes C Centrifuge tubes, clear Isotonic saline Pipettes Ceph floc antigen Cotton Lab record book	A Icteus index 1. Dilute serum 2. Read O.D. of standard and test 3. Calculate value B Thymol turbidity 1. Mix reagents 2. Read turbidity C Ceph Floc 1. Mix reagents 2. Let stand in dark 3. Read flocculation	Safety Reliability of results Chemical handling Glassware handling Sources of error Hazard Patient's well-being Contamination Lacerations
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

(TASK STATEMENT) PERFORM ICTEUES INDEX, THYMOL TURBIDITY, CEPHALIN CHOLESTEROL FLOCCULATION

SCIENCE	MATH - NUMBER SYSTEMS
Digestion and metabolism of foods Liver function and diseases Red cell breakdown Chemical reactions of tests Normal values Interpretation of results Bile formation and function	OD unknown OD standard x 10 = Icteuues index Dilutions Basic arithmetic skills Ratios and proportions Measure of time Metric measure of volume
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
Reading Writing	Procedure Report results, label, draw graph
	SKILLS/CONCEPTS Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage
	161

(TASK STATEMENT) PERFORM BSP (BROMSULFONPHTHALIN) TEST

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>BSP dye Serum sample Venipuncture equipment Sodium hydroxide Hydrochloric acid Pipettes Spectrophotometer Cuvettes Distilled water Graph paper Lab record book Timer</p>	<p>Calculate volume dye needed (mg/kg body weight) Administer dye Draw blood after 45 - 60 seconds Develop color Calibrate spectrophotometer Read % T Determine mg % dye from graph</p>	<p>Safety Reliability and accuracy of results Proper glassware handling Chemical handling Proper amount dye Hazards Patient's well-being Lacerations</p>
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> Technician does not inject dye generally Timing accuracy is essential Blood sample should be drawn from opposite arm from where dye injected</p>	<p><u>ERRORS</u> Incorrect results</p>

(TASK STATEMENT) PERFORM BSP (BROMSULFONPHTHALIN) TEST

SCIENCE	MATH - NUMBER SYSTEMS
<p>Function of BSP dye Liver function Red cell breakdown Liver diseases Bile formation and function Digestion and metabolism Chemical reactions of tests</p>	<p>Basic arithmetic skills Measure with the Metric and English system and convert between them Ratio and proportions [mg dye/kilogram body weight] Measures of weight Measure of metric volume Measure of time Read graph</p>
COMMUNICATIONS	
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>
<p>Reading Writing</p>	<p>Procedure Report results, label, draw graph</p> <p>163</p>
<u>SKILLS/CONCEPTS</u>	
<p>Comprehension, medical terminology, description of mechanism Penmanship, spelling, accuracy, progress report, usage</p>	

Duty G Performing Urinalysis

- 1 Collect and preserve urine specimens
- 2 Determine physical characteristics of urine
- 3 Determine pH of urine
- 4 Determine specific gravity of urine
- 5 Determine glucose in urine
- 6 Determine protein in urine
- 7 Determine presence of ketone bodies in urine
- 8 Perform test for bile and urobilinogen in urine
- 9 Perform test for blood in urine and feces
- 10 Perform microscopic examination of urine
- 11 Perform renal function test

170

(TASK STATEMENT) COLLECT AND PRESERVE URINE SPECIMENS

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Clean urine bottle or cup and lid	Collect random specimen (voided or catheterized) Collect 24 hour specimen Preserve 24 hour specimen refrigeration toluene thymol formalin	Hazard Contamination
<u>DECISIONS</u> Determine standards for accuracy	<u>CUES</u> Test requested determines time specimen taken Fresh specimen use prevents decomposition Unless specimen is preserved, test should be run within 1-2 hours Catheterized specimens needed for bacteriologic examination First morning specimen is usually best	<u>ERRORS</u> Poor quality sample Contamination

(TASK STATEMENT) COLLECT AND PRESERVE URINE SPECIMENS

SCIENCE		MATH -- NUMBER SYSTEMS
Composition of urine Anatomy and physiology of urinary system Normal quantity of urine Effects of diet on urine		
COMMUNICATIONS		
<u>PERFORMANCE MODES</u> Reading	<u>EXAMPLES</u> Requisitions	<u>SKILLS/CONCEPTS</u> Comprehension, medical terminology
	167	172

(TASK STATEMENT) DETERMINE PHYSICAL CHARACTERISTICS OF URINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Urine specimen Laboratory record book Graduated glassware</p>	<p>Determine volume Observe and distinguish color Assess transparency Determine odor</p>	<p>Safety Wash hands Hazard Contamination</p>
<p><u>DECISIONS</u></p> <p>Determine need for further testing</p>	<p><u>CUES</u></p> <p>Physical characteristics can give clues to findings in other portions of urinalysis When abnormal characteristics are found, certain other tests are indicated</p>	<p><u>ERRORS</u></p> <p>Unnecessary tests Fail to perform needed tests</p>

(TASK STATEMENT) DETERMINE PHYSICAL CHARACTERISTICS OF URINE

SCIENCE	MATH — NUMBER SYSTEMS
<p>Anatomy and physiology of urinary system Causes and characteristics of abnormal results Quantity of urine excretion Urinary pigments</p>	<p>Metric system measure of volume</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u> Reading Viewing</p>	<p><u>EXAMPLES</u> Results Specimen</p> <p>169</p>
<p><u>SKILLS/CONCEPTS</u> Comprehension, technical terminology Visual analysis, Color discrimination</p>	

(TASK STATEMENT) DETERMINE pH OF URINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Urine specimen Lab record book Combistix Bililabstix Nitiazine paper Reference color charts</p>	<p>Use combistix (ames) Use Bililabstix (ames) Use nitrazine paper Compare with color charts</p>	
<p><u>DECISIONS</u></p> <p>Determine standards for accuracy</p>	<p><u>CUES</u></p> <p>pH changes after specimen stands pH can give clue to microscopic identification Alkaline pH can cause dissolving of some elements</p>	<p><u>ERRORS</u></p> <p>Improper results</p>

(TASK STATEMENT) DETERMINE pH OF URINE

SCIENCE	MATH — NUMBER SYSTEMS
<p>Kidney function Hydrogen ion concentration Acid-base theories Effects of pH on body Normal = 5-7</p>	<p>Scientific notation</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing Viewing</p>	<p><u>EXAMPLES</u></p> <p>Results Color charts</p> <p>171</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology Penmanship, detail/inference Visual analysis, Color discrimination</p>

(TASK STATEMENT) DETERMINE SPECIFIC GRAVITY OF URINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Urine specimen Urinometer Hydrometer jar Lab record book Clinical refractometer Distilled water Thermometer</p>	<p>Calibrate urinometer Correct for temperature Read specific gravity on urinometer Correct for glucose present Use refractometer</p>	<p>Safety Urinometers are fragile Contain mercury Hazard Mercury vapors Lacerations</p>
<p><u>DECISIONS</u></p> <p>Select urinometer Determine standards for accuracy</p>	<p><u>CUES</u></p> <p>Bubbles in urine can cause errors in reading Urinometer must float free off the bottom of the container Read on flat surface at eye level Read at bottom of meniscus</p>	<p><u>ERRORS</u></p> <p>No reading</p>

(TASK STATEMENT) DETERMINE SPECIFIC GRAVITY OF URINE

SCIENCE	MATH - NUMBER SYSTEMS	
Kidney function Density Specific gravity Formation of solutions Urinary physiology Relationship of volume and pH Osmolality Refractive index Meniscus	$\frac{\text{Weight of solution}}{\text{Weight of water}} = \text{Specific gravity}$ Temperature correction = .001/degree above or below calibration level Glucose correction = .003/1 gm glucose Temperature measurement	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing Viewing	Report results, note correction Specimen	Comprehension, medical terminology Progress report, spelling, terminology, logic Visual analysis, Detail/inference

(TASK STATEMENT) DETERMINE GLUCOSE IN URINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Urine specimen Clinitest tablets Tes-tape Clinistix Bililabstix Benedict's qualitative reagent (copper sulfate, sodium citrate, sodium carbonate, distilled water) Test tube Boiling water bath Reference color charts Lab record book Cold water bath Timer Water pipets Quantitative test: Benedict's quantitative reagent Sodium carbonate Heat source Pipet Test tubes	Non-specific tests - reducing activity Benedict's qualitative test Clinitest (Ames) Specific tests - using enzyme oxidase Tes-tape (Eli Lilly) Clinistix (Ames) Bililabstix (Ames) or Combistix (Ames) Report intensity of color change Quantitative tests (24 hour specimen) Benedict's quantitative Read results	Safety Keep bottles of stix and tablets tightly closed Hazard Deterioration of chemicals
Determine standards for accuracy	Accurate timing is critical Read procedures carefully Clinitest can pass through 4 plus results to show higher value; watch as reaction occurs	Incorrect results
<u>DECISIONS</u>	<u>CUES</u>	<u>ERRORS</u>

(TASK STATEMENT) DETERMINE GLUCOSE IN URINE

SCIENCE	MATH - NUMBER SYSTEMS
<p>Carbohydrate digestion and assimilation Chemical reactions in testing Diabetes - causes, symptoms, diagnosis, results Kidney function - absorption Renal threshold concept</p>	<p><u>Quantitative glucose</u></p> $\frac{1}{\text{cc urine}} = \% \text{ glucose}$ $\frac{\text{Total volume}}{\text{cc urine}} = \text{gms/24 hours}$ <p>100</p> <p>Metric system measure of volume Basic arithmetic skills</p>
COMMUNICATIONS	
PERFORMANCE MODES	EXAMPLES
<p>Reading Writing</p>	<p>Procedures supplied with tablets and sticks carefully Record results</p> <p>175</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, detail/inference, medical terminology Penmanship, spelling, logic</p>

(TASK STATEMENT) DETERMINE PROTEIN IN URINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Urine specimen Centrifuge Centrifuge tubes Water Sodium chloride Glacial acetic acid Boiling water bath Timer Cold water bath Bumintest tablets Sodium sulfate Sulfsalicylic acid Standards Test tubes Albutest tablets Albustix Combistix Billabstix</p>	<p>Perform heat and acetic acid test Perform Bumintest (Ames) Perform Sulfsalicylic acid test Perform Albutest (Ames) Perform Albustix (Ames) Perform Combistix or Billabstix Read results</p>	
<p><u>DECISIONS</u></p> <p>Turbidity comparison</p>	<p><u>CUES</u></p> <p>Presence of protein can indicate findings in sediment Mucin interferes with protein tests</p>	<p><u>ERRORS</u></p>

(TASK STATEMENT) DETERMINE PROTEIN IN URINE

(TASK STATEMENT) DEPENDENT PROBLEM IN ORIGIN		MATH - NUMBER SYSTEMS	
SCIENCE			
Kidney function Filtration and reabsorption - urine formation Diseases causing protein occurrence Chemical reactions of protein tests Protein digestion and assimilation			
COMMUNICATIONS			
<u>PERFORMANCE MODES</u> Reading Writing		<u>EXAMPLES</u> Literature supplied by manufacturer Record results	<u>SKILLS/CONCEPTS</u> Comprehension, medical terminology Penmanship, spelling
		177	

(TASK STATEMENT) DETERMINE PRESENCE OF KETONE BODIES IN URINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Urine specimen Sodium nitroprusside Ammonium sulfate Ammonium hydroxide Test tubes Acetest tablets Ketostix 10% ferric chloride Timer Heat source</p>	<p>Perform Rothera's test Use acetest (Ames) Use Ketostix (Ames) Perform Jerhardt's test Read results</p>	
<p><u>DECISIONS</u> Determine standards for accuracy</p>	<p><u>CUES</u> These tests are performed when glucose is positive Timing is critical</p>	<p><u>ERRORS</u> Incorrect results</p>

(TASK STATEMENT) DETERMINE PRESENCE OF KETONE BODIES IN URINE

SCIENCE	MATH - NUMBER SYSTEMS
<p>Digestion and assimilation of fats Diabetes mellitus causes and effects Types of ketone bodies Physiologic effects of ketone accumulation Chemical reactions of tests</p>	<p>Dilutions Metric system measure of volume</p>
COMMUNICATIONS	
PERFORMANCE MODES	SKILLS/CONCEPTS
<p>Reading Writing</p>	<p>Comprehension, medical terminology Penmanship, spelling</p>
<p>Literature supplied by manufacturer Record results</p>	<p>179</p>

(TASK STATEMENT) PERFORM TEST FOR BILE AND UROBILINOGEN IN URINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Urine specimen Lab record book Fouchet's reagent Barium chloride paper Ictotest (Ames) Water Color comparison chart Ehrlich's reagent Sodium acetate Chloroform Butanol Test tube Separatory funnel Graduated glassware Billilabstix</p>	<p>Perform Harrison's test Use Ictatest tablets (Ames) Perform Ehrlich's qualitative test for urobilinogen Use Billilabstix (Ames) Read results</p>	<p>Safety Chemicals used in Ehrlich's test chlorinated are aromatic solvents Hazard Inhalation of chemical fumes</p>
<p><u>DECISIONS</u></p> <p>Select timing device</p>	<p><u>CUES</u></p> <p>Color of urine indicates need to perform test for bile Timing is critical</p>	<p><u>ERRORS</u></p> <p>Misread results Improper timing</p>

(TASK STATEMENT) PERFORM TEST FOR BILE AND UROBILINOGEN IN URINE

SCIENCE	MATH - NUMBER SYSTEMS
<p>Kidney function Digestive system RBC breakdown process Chemical reactions of RBC destruction Chemical reactions of tests for bile Liver function</p>	<p>Metric system measure of volume</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Literature supplied by manufacturer Record results</p> <p>181</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology Penmanship, spelling</p>

(TASK STATEMENT) PERFORM TEST FOR BLOOD IN URINE AND FECES

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Urine specimen or feces Lab record book Glacial acetic acid Hydrogen peroxide Gum guaiac Ethanol Filter paper Dropping tube Hemastix Hematest Color comparison charts Bililabstix</p>	<p>Perform gum guaiac test Use Hemastix (Ames) Use Her test tablets (Ames) Use Bililabstix (Ames) Read results</p>	<p>Safety Use chemicals with care Hazard Burr.s</p>
<p><u>DECISIONS</u></p> Identify blood	<p><u>CUES</u></p> Many chemical substances, if present, interfere with these tests	<p><u>ERRORS</u></p> Misread sample

(TASK STATEMENT) PERFORM TEST FOR BLOOD IN URINE AND FECES

SCIENCE	MATH — NUMBER SYSTEMS
Kidney function RBC morphology and breakdown Chemical reactions of tests	Metric measure of volume
COMMUNICATIONS	
<u>PERFORMANCE MODES</u> Reading Writing	<u>SKILLS/CONCEPTS</u> Comprehension, medical terminology Penmanship, spelling
<u>EXAMPLES</u> Literature supplied by manufacturer Record results	183

(TASK STATEMENT) PERFORM MICROSCOPIC EXAMINATION OF URINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Urine specimen Centrifuge Centrifuge tubes Glass slides Cover slip Microscope Lab record book</p>	<p>Centrifuge specimen Decant supernatant Place sediment on slide Examine microscopically Identify normal and abnormal constituents under high and lower power and estimate number</p>	<p>Hazard Scratched or cracked lens Broken slides</p>
<p><u>DECISIONS</u></p> <p>Select equipment and procedure Identify type of crystals</p>	<p><u>CUES</u></p> <p>Remove amorphous before centrifugation Standing of urine cause distortion or breakdown of sediment</p>	<p><u>ERRORS</u></p> <p>Misread sample</p>

(TASK STATEMENT) PERFORM MICROSCOPIC EXAMINATION OF URINE

SCIENCE	MATH - NUMBER SYSTEMS
<p>Kidney physiology and anatomy Digestion process Characteristics, morphology of biologic and chemical sediment Blood cells Casts Crystals Epithelial cells Bacteria Miscellaneous substances</p>	<p>Counting sequentially</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing Viewing</p>	<p><u>EXAMPLES</u></p> <p>Record results Specimen</p> <p>185</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology Penmanship, spelling, accuracy Visual analysis, Detail/inference</p>

(TASK STATEMENT) PERFORM RENAL FUNCTION TEST

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Properly collected urine specimens 1000 ml graduated cylinder Water Sodium hydroxide Spectrophotometer Graph paper Urinometer Hydrometer jar Lab record book</p>	<p>Perform PSP (Phenolsulf on phthalein) test Measure volume Adjust pH Calibrate spectrophotometer Read % T of unknowns Determine value from graph Perform Mosenthal's Concentration Test Measure volume Measure specific gravity</p>	<p>Safety Sodium hydroxide is corrosive Hazard Burns</p>
<p><u>DECISIONS</u></p>	<p><u>CUES</u></p>	<p><u>ERRORS</u></p>

(TASK STATEMENT) PERFORM RENAL FUNCTION TESTS

SCIENCE	MATH - NUMBER SYSTEMS
<p>Kidney function Absorption and filtration Liver function Hydrogen ion concentration Water balance in body Specific gravity and density Kidney pathology</p>	<p>Draw graph Metric measure of volume Basic arithmetic skills</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Graph Record results</p> <p>187</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology Penmanship, spelling, accuracy</p> <p>192</p>

Duty H Performing Bacteriology Procedures

- 1 Sterilize equipment and supplies
- 2 Prepare culture media
- 3 Inoculate culture media
- 4 Incubate cultures
- 5 Determine culture morphology
- 6 Prepare bacterial slides
- 7 Stain slides using gram stain
- 8 Determine microscopic morphology
- 9 Identify bacteria
- 10 Determine drug sensitivity

(TASK STATEMENT) STERILIZE EQUIPMENT AND SUPPLIES

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY — HAZARD
<p>Equipment Glassware Used or contaminated equipment Autoclave Burner (Bunsen or alcohol) Water source 75% alcohol Soap solution Internal timer</p>	<p>Sterilize transfer needles and loops by burning Wrap small equipment and supplies Autoclave at 121°C for 15 - 20 minutes at 15 pounds pressure</p>	<p>Safety Open fire Pressure in autoclave and heat Glassware can be broken Improper sterilization will not kill bacteria Hazard Burns Explosion Lacerations Bacterial contamination</p>
DECISIONS	CUES	ERRORS
<p>Select methods of sterilization Select autoclave</p>	<p>Method used depends on nature of material Use of autoclave (steam under pressure) is most effective</p>	<p>Contamination</p>

(TASK STATEMENT) . STERILIZE EQUIPMENT AND SUPPLIES

SCIENCE	MATH — NUMBER SYSTEMS
<p>Resistance of microorganisms to sterilization Gas laws Effect of chemicals on microorganisms</p>	<p>Measure of temperature [to include Kelvin] Ratio and proportion</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Thermometer, gauge</p> <p>191</p>
<p><u>SKILLS/CONCEPTS</u></p> <p>Detail/inference</p>	

TASK STATEMENT) PREPARE CULTURE MEDIA

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY — HAZARD	196
<p>Dehydrated media Distilled water Balance Weighing paper or dish Volumetric glassware Autoclave Water source Petri dishes Sterile tubes Plugs Nutrient additives, i.e. blood Heat source Beakers Stirring device</p>	<p>Accurately weigh correct amount of dehydrated media Dilute accurately Sterilize Pour plates or tubes agar slants broth tubes agar plates fermentation tubes Use sterile transfer technique Store in refrigerator</p>	<p>Safety Use sterile techniques Hazard Burn Lacerations Contaminated plates</p>	
<p><u>DECISIONS</u></p> <p>Determine type media necessary</p>	<p><u>CUES</u></p> <p>Use of media Shelf-life of prepared media Accurate timing is essential Denaturation of sugars at elevated temperature and pressure</p>	<p><u>ERRORS</u></p> <p>Media will not set-up</p>	

(TASK STATEMENT) PREPARE CULTURE MEDIA

SCIENCE	MATH - NUMBER SYSTEMS
<p>Cultural requirements of bacteria Selective, differential and enrichment media characteristics</p>	<p>Proportions and ratios Percent measures Measures of metric volume Measure of metric weight</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u> Reading</p>	<p><u>EXAMPLES</u> Directions for preparation on media</p> <p><u>SKILLS/CONCEPTS</u> Comprehension, medical terminology</p> <p>193</p>

(TASK STATEMENT) INOCULATE CULTURE MEDIA

198

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Culture specimen Media Bunsen burner or alcohol lamp Transfer loop or needle	Determine appropriate culture media Sterilize loop or needle, tube lips Use sterile technique for transfer Use appropriate inoculation techniques	Safety Sterile techniques must be utilized Hazard Burns Bacterial contamination
<u>DECISIONS</u> Select inoculation technique	<u>CUES</u> Media used dependent on type specimens Inoculation techniques is dependent on type media	<u>ERRORS</u> Improper culture growth

(TASK STATEMENT) INNOCULATE CULTURE MEDIA

SCIENCE	MATH — NUMBER SYSTEMS
<p>Normal and abnormal flora Characteristics of selective, differential, and enrichment media Cultural requirements of bacteria</p>	
<p>COMMUNICATIONS</p>	
<p><u>PERFORMANCE MODES</u> Reading</p>	<p><u>EXAMPLES</u> Difco manual</p>
	<p><u>SKILLS/CONCEPTS</u> Comprehension, terminology</p>
	<p>195</p>
	<p>199</p>

(TASK STATEMENT) INCUBATE CULTURES

<p>TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON</p>	<p>PERFORMANCE KNOWLEDGE</p>	<p>SAFETY -- HAZARD</p>
<p>Properly inoculate culture media Incubator Brerwer anaerobic jar</p>	<p>Incubate under aerobic conditions Incubate using anaerobic conditions Brerwer anaerobic jar Innoculate deep in solid media Allow time for growth</p>	<p>Hazard Electrical shocks</p>
<p><u>DECISIONS</u> Select time and temperature</p>	<p><u>QUES</u> Proper temperature must be controlled Incubator must contain moisture to avoid drying out of cultures</p>	<p><u>ERRORS</u> Ruin culture</p>

200

(TASK STATEMENT) INCUBATE CULTURES

SCIENCE	MATH -- NUMBER SYSTEMS
<p>Cultural requirements of bacteria Aerobic and anaerobic conditions Optimum growth requirements</p>	<p>Temperature measurement</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Difco manual</p>
<p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, terminology</p>	<p>197</p>

(TASK STATEMENT) DETERMINE CULTURE MORPHOLOGY

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Mature culture Lab record book Hand lens or dissecting microscope Calory counter</p>	<p>Observe color Observe characteristics of colony, i.e. texture, shape, luminous characteristics Observe hemolysis on blood agar Observe motility Observe odor Count calories</p>	<p>Safety Avoid contamination Hazard Contamination of lab and culture</p>
<p><u>DECISIONS</u></p>	<p><u>CUES</u></p>	<p><u>ERRORS</u></p>

(TASK STATEMENT) DETERMINE CULTURE MORPHOLOGY

SCIENCE	MATH — NUMBER SYSTEMS
Cultural characteristics of bacteria	
COMMUNICATIONS	
<u>PERFORMANCE MODES</u> Reading Writing Viewing	<u>EXAMPLES</u> Record results Mature culture characteristics
	<u>SKILLS/CONCEPTS</u> Comprehension, terminology Penmanship, spelling Visual analysis, Detail/inference, Color discrimination
	199
	203

(TASK STATEMENT) PREPARE BACTERIAL SLIDES

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Slides Mature culture Wire loop or needle Bunsen burner or alcohol burner Distilled water</p>	<p>Make direct smears from specimen Make smears from both cultures Emulsify growth from plate in water on slide Air-dry and fix with heat</p>	<p>Safety Use proper sterile technique Hazards Burns Lacerations from slides Contamination of lab and culture</p>
<p><u>DECISIONS</u></p> <p>Select procedure and storage</p>	<p><u>CUES</u></p> <p>Excessive heating can destroy bacteria Single layer sample desirable</p>	<p><u>ERRORS</u></p> <p>Lost slide</p>

TASK STATEMENT) PERPARE BACTERIAL SLIDES

TASK STATEMENT: PREPARE BACTERIAL SLIDES		MATH — NUMBER SYSTEMS	
SCIENCE			
Effects of heat on microorganisms			
COMMUNICATIONS			
<u>PERFORMANCE MODES</u> Reading Writing	<u>EXAMPLES</u> Labels	<u>SKILLS/CONCEPTS</u> Comprehension, terminology Penmanship, spelling, accuracy	205
	201		

(TASK STATEMENT) STAIN SLIDES USING GRAM STAIN

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Fixed slide Staining rack Timer Blotting paper Crystal violet Ethyl alcohol Ammonium oxalate Distilled water Iodine solution Safranin Acetone Squeeze bottles Microscope Immersion oil Lab record book	Follow gram-stain procedures Air dry Examine microscopically under oil immersion	Safety Stains are toxic and flammable Wear protective devices Hazard Fumes Dyes stain hand and clothes Fire
<u>DECISIONS</u> Determine standards for accuracy	<u>CUES</u> Proper timing is essential Thorough washing necessary Avoid direct contact with slides Culture must be within the 18 - 24 hour period	<u>ERRORS</u> Poor quality slide Over stained Under stained

TASK STATEMENT) STAIN SLIDES USING GRAM STAIN

SCIENCE	MATH — NUMBER SYSTEMS	
Gram positive and gram negative differentiation Chemical reactions of stains with bacteria	Measurement of time Basic arithmetic skills	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Staining procedure Record results	Comprehension, terminology Penmanship, spelling, accuracy
	203	207

(TASK STATEMENT) DETERMINE MICROSCOPIC MORPHOLOGY

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Stained slide Microscope Immersion oil Lab record book</p>	<p>Classify gram reaction (gram + or -) Observe size shape arrangement of cells internal cellular structure</p>	<p>Safety Use microscope correctly Hazard Cracked or scratched microscope lens</p>
<p><u>DECISIONS</u> Identify proper sample</p>	<p><u>CUES</u> Visual observation</p>	<p><u>ERRORS</u> Misread sample</p>

209

(TASK STATEMENT) IDENTIFY BACTERIA

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY – HAZARD 210
<p>Selective and differential media Culture specimen Culture characteristics Microscopic characteristics</p>	<p>Determine culture morphology Determine microscopic morphology Determine staining characteristics Determine chemical reactivity Determine differential media reaction Determine optimum culturing temperature Organize flow chart Classify</p>	<p>Safety Use sterile technique Hazard Contamination of lab and culture</p>
<p><u>DECISIONS</u> Classify bacteria</p>	<p><u>CUES</u> Bacteria patterns</p>	<p><u>ERRORS</u> Misread sample</p>

ASK STATEMENT) IDENTIFY BACTERIA

SCIENCE		MATH — NUMBER SYSTEMS	
Classification of bacteria Characteristics of bacteria Normal flora Pathological organisms			
COMMUNICATIONS			
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS	
Reading Writing Viewing	Difco manual Handbook of microbiology Record results Characteristics of bacteria	Comprehension, technical terminology, vocabulary Penmanship, spelling, accuracy Visual analysis, Detail/inference	
	207	211	

(TASK STATEMENT) DETERMINE DRUG SENSITIVITY

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY — HAZARD
Pure culture Sensitivity discs	Streak plate with pure culture Add sensitivity discs Incubate Read zone of sensitivity to drugs	Safety Use sterile techniques Hazard Contamination of culture and lab
<u>DECISIONS</u> Determine type of sensitivity disc	<u>CLUES</u> Type of drug	<u>ERRORS</u> Improper results

ASK STATEMENT) DETERMINE DRUG SENSITIVITY

SCIENCE	MATH -- NUMBER SYSTEMS	
Bacterial reaction to drugs Normal flora	Use metric measure of length	
COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Reading Writing	Reference material Record results	Comprehension, technical terminology Penmanship, spelling
	209	

Duty I Performing Blood Bank and Serology Procedures

- 1 Perform ABO typing**
- 2 Perform Rh typing**
- 3 Crossmatch**
- 4 Perform VDRL**

214

(TASK STATEMENT) PERFORM ABO TYPING

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Saline - isotonic Blood sample Anti-A typing sera Anti-B typing sera Slide Test tubes Microscope Centrifuge Lab record book Pipette A cells B cells</p>	<p>Place blood samples in tray Add antigen Observe results</p>	<p>Safety Typing of blood is a life or death procedure Hazard Transfusion reaction Death to patient</p>
<p><u>DECISIONS</u> Determine type of test technique</p>	<p><u>CUES</u> Type of test</p>	<p><u>ERRORS</u> Contaminated or old typing sera Incorrect cell suspension Timing insufficient Identification insufficient Low sensitivity Improper reading of agglutination Contaminated specimens Glassware dirty</p>

TASK STATEMENT) PERFORM ABO TYPING

SCIENCE		MATH - NUMBER SYSTEMS
<p>Volume of blood Composition of blood History of transfusions Anticoagulant theory Inheritance of blood groups Antigen - antibody reactions Blood group systems Means of detecting antigen - antibody reactions ABO blood group system Titer Agglutination Typing sera characteristics</p>		Proportions
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing Viewing	Reference manuals, procedures Record results, labels Observe results 213	Comprehension, medical terminology Penmanship, spelling, accuracy Visual analysis, Detail/inference 246

(TASK STATEMENT) PERFORM Rh TYPING

217

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Anti-Rh typing Blood sample - whole blood Test tubes Isotonic saline Incubator or water bath Microscope Pipette Slide Warning box Lab record book Timer	Slide type Tube type	Safety Typing of blood is a life or death procedure Hazard Transfusion reaction Death to patient
<u>DECISIONS</u> Determine type of test technique	<u>CUES</u> Type of test	<u>ERRORS</u> Contaminated or old typing serum Incorrect temperature Incorrect cell suspension Confusion of rouleaux formation with agglutination Contaminated specimen Insufficient identification Glassware dirty

217

TASK STATEMENT) PERFORM Rh TYPING

SCIENCE	MATH - NUMBER SYSTEMS
<p>Composition of blood History of transfusions Anticoagulant theory Inheritance of blood groups Antigen - antibody reactions Rh blood group system and theory Agglutination Erythroblastosis fetalis</p>	<p>Measurement of temperature Proportions</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading Writing</p>	<p><u>EXAMPLES</u></p> <p>Reference manuals, procedures Record results, labels</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, medical terminology Penmanship, spelling, accuracy</p>

(TASK STATEMENT) CROSSMATCH

249

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
Patient serum and cells Donor serum and cells Coomb's serum Test tubes Pipettes Incubator and waterbath Centrifuge Slide Microscope Lab record book Wash bottle Saline Albumin	Perform ABO typing Perform Rh typing Select matching unit Set up crossmatch saline high-protein coombs	Safety Clerical errors Organizational errors Technical errors Hazard Transfusion reactions Patient death
<u>DECISIONS</u> Determine standards for accuracy	<u>CUES</u> Detects - irregular antibodies incompatibilities labeling, identifying errors Work should be checked by another technician Timing is critical	<u>ERRORS</u> Incorrect matching

ASK STATEMENT) CROSSMATCH

SCIENCE		MATH - NUMBER SYSTEMS
Crossmatch theory Antigen - antibody reactions Blood group systems Composition of blood Agglutination Causes of error Glassware cleaning Selection of donors Blood derivatives Clerical work in blood bank Transfusion reactions Coomb's reaction		Ratio - proportion Measure of temperature Measure of time
COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Reading Writing	Reference manuals, procedures Record results, labels 217	Comprehension, medical terminology Penmanship, spelling, accuracy

(TASK STATEMENT) PERFORM VDR

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Serum samples Controls Pipettes Ring slides Flat bottom, glass stoppered bottle Syringe and needle Rotator Waterbath or incubator Microscope Lab record book Buffered saline Antigen Timer</p>	<p>Prepare antigen Inactivate serum Pipette serum and antigen Rotate Examine for agglutination Run positive and negative controls</p>	<p>Safety Use serum cautiously Hazard Technician contamination</p>
<p><u>DECISIONS</u></p> <p>Select rotation time, speed and temperature Select antigen</p>	<p><u>CUES</u></p> <p>Rotation time and speed must be exact Antigen must be prepared properly Accurate measuring essential</p>	<p><u>ERRORS</u></p> <p>Misread sample Wrong antigen</p>

SCIENCE		MATH - NUMBER SYSTEMS
Antigen - antibody reactions Tests for syphilis Titers Immunity - resistance to disease Mechanism of reaction False negative and positive results Pathological aspects of syphilis		Ratios and proportions Measurement of volume
COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Reading Writing Viewing	Reference manuals, procedures Record results, labels Serum for agglutination	Comprehension, medical terminology Penmanship, spelling, accuracy Visual analysis, Detail/inference
	219	222

Duty J Performing Histology Procedures

- 1 Fix tissues**
- 2 Embedded tissue in paraffin**
- 3 Cut thin sections and mount on slides**
- 4 Stain tissue slides**
- 5 Prepare frozen sections**

223

(TASK STATEMENT) FIX TISSUES

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY — HÁZARD 224
Tissue specimen Auto technicon Fixing solutions	Place appropriate size tissue in auto technicon Set time limits for each solution and temperature	Safety Fixing agents must be handled correctly Hazard Inhalation of fumes Chemical burns
<u>DECISIONS</u> Select fixation agent Determine standards for accuracy	<u>CUES</u> Nature of tissue, stain to be used Process must be exact to obtain proper penetration, preserving and harden- ing of tissues Tissue blocks must be fairly small to allow proper fixation	<u>ERRORS</u> Improper penetration

SCIENCE	MATH - NUMBER SYSTEMS
<p>Tissue preservation theory Properties of fixative agents</p>	<p>Measurement of time</p>
COMMUNICATIONS	
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Operating manual</p> <p>223</p>
<p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, technical terminology</p> <p>225</p>	

(TASK STATEMENT) EMBEDDED TISSUE IN PARAFFIN

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Fixed tissue sample Paraffin Heat source Pan for embedding Cold water Forceps Knife Pencil</p>	<p>Melt paraffin Transfer tissue to paraffin Cool rapidly Trim blocks Mark blocks</p>	<p>Hazard Paraffin burns</p>
<p><u>DECISIONS</u></p> <p>Determine type of media</p>	<p><u>CUES</u></p> <p>Paraffin melts at 56° White patches are caused by clearing agent carried over from fixing process</p>	<p><u>ERRORS</u></p> <p>Ruin sample</p>

SCIENCE		MATH - NUMBER SYSTEMS
COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Reading Writing	Procedures Labels	Comprehension, terminology Penmanship, spelling, accuracy
	225	227

(TASK STATEMENT) CUT THIN SECTIONS AND MOUNT ON SLIDES

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Embedded tissue Microtome Warm water bath Slides Fixative Glass pick</p>	<p>Cut sections on microtome Float sections in warm water Mount on slide with fixative Drain and dry Label</p>	<p>Safety Microtome knives are extremely sharp Hazards Lacerations Destruction of tissue</p>
<p><u>DECISIONS</u></p> <p>Determine standards for accuracy</p>	<p><u>CUES</u></p> <p>Errors are usually due to one or more causes: paraffin temperature, sharpness of knife, angle of knife</p>	<p><u>ERRORS</u></p> <p>Poor mount Unusable slide</p>

(TASK STATEMENT) CUT THIN SECTIONS AND MOUNT ON SLIDES

SCIENCE		MATH — NUMBER SYSTEMS	
Surface tension Capillary adhesion Errors of cutting sections		Metric measurement of length	
COMMUNICATIONS			
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>	
Reading Writing	Label slide	Comprehension, medical terminology Penmanship, spelling, accuracy	
	227	229	

(TASK STATEMENT) STAIN TISSUE SLIDES

230

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Properly mounted slide Staining trays Staining jars Hematoxylin Xylol Alcohol Distilled water HCL Lithium carbonate Eosin Carbo-xylol Labels Cover slips	Follow staining procedure Mount with cover slip Label completely	Safety Use stains carefully Use glassware properly Hazards Stains dye clothes and skin Chemical inhalation Lacerations
<u>DECISIONS</u> Determine type of stain and time	<u>CUES</u> Accurate timing is essential Tissue characteristics	<u>ERRORS</u> Ruin sample Poor quality slide

TASK STATEMENT) STAIN TISSUE SLIDES

SCIENCE		MATH — NUMBER SYSTEMS
Color index for dyes Chemical reactions of staining Affinity of cells parts to stain Regressive and progressive staining	Measurement of time	
COMMUNICATIONS		
<u>PERFORMANCE MODES</u>	<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Reading Writing Viewing	Procedure Label slide Stained slide 229	Comprehension, medical terminology Penmanship, spelling, accuracy Visual analysis, Color discrimination

(TASK STATEMENT) PREPARE FROZEN SECTIONS

232

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Tissue Heat source Formalin Freezing microtome Water Carbon monoxide Pan of water Slide Glass pick Alcohol Hamatoxylin Eosin Xylol Coverslip Label	Boil tissue in formalin Freeze Cut section Float onto slide Seal with heat Stain Coverslip Mount	Safety Use heat carefully around volatile solutions Use care in cutting Use glassware correctly Hazard Burns Fire Lacerations
<u>DECISIONS</u> Determine quality of sample	<u>CUES</u> Done for rapid preparation Done for special stains	<u>ERRORS</u> Section too thick

ASK STATEMENT) PREPARE FROZEN SECTIONS

SCIENCE	MATH — NUMBER SYSTEMS	
Fixation process Staining reactions Freezing process - chemical reactions	Measurement of time	
COMMUNICATIONS		
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS
Reading Writing	Procedure Label slide	Comprehension, medical terminology Penmanship, spelling, accuracy

2008

201

Duty K Performing EKG

- 1 Prepare patient for EKG**
- 2 Set up electrocardiograph**
- 3 Operate electrocardiograph**
- 4 Perform follow-up care of patient**
- 5 Perform follow-up care of machines**
- 6 Mount electrocardiogram**

(TASK STATEMENT) PREPARE PATIENT FOR EKG

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Patient Warm, quiet room Bed Electrode jelly Strap electrodes Suction electrode</p>	<p>Reassure patient Situuate patient in comfortable position Apply electrode jelly Fasten electrode straps to limbs Apply suction electrode to first chest position</p>	<p>Safety Patient should not be touching any metal objects</p>
<p><u>DECISIONS</u></p> <p>Determine proper hook-up Locate lead hook-up</p>	<p><u>CUES</u></p> <p>Patient must lie still and avoid talking Excess jelly extorts results Attach electrode to stump in case of amputation Straps should provide firm contact, but not discomfort</p>	<p><u>ERRORS</u></p> <p>Improper reading</p>

ASK STATEMENT) PREPARE PATIENT FOR EKG

SCIENCE		MATH — NUMBER SYSTEMS	
Composition of electrode jelly Purpose of EKG			
COMMUNICATIONS			
PERFORMANCE MODES	EXAMPLES	SKILLS/CONCEPTS	
Reading Speaking	Requisition slip Reassure patient	Comprehension Appropriate diction, clarity of expression, persuasion, gestures, poise	

(TASK STATEMENT) SET UP ELECTROCARDIOGRAPH

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY — HAZARD
<p>Prepared patient Electrocardiograph Cable cards Paper for tracing</p>	<p>Ground machine Allow adequate warm-up Attach cable cards to patient electrode straps Check paper</p>	<p>Safety Ground machine Attach leads correctly Hazard Electrical shock Interference</p>
<p><u>DECISIONS</u></p> <p>Determine lead hook-up Proper operation of machine</p>	<p><u>CUES</u></p> <p>Mix-up in lead attachments causes abnormalities in tracing Proper grounding eliminates A-C interference Oxygen tents may need to be turned off by nurse</p>	<p><u>ERRORS</u></p> <p>Inaccurate reading</p>

ISK STATEMENT) SET UP ELECTROCARDIOGRAPH

SCIENCE	MATH — NUMBER SYSTEMS
Theory of machine operation -electronic amplifier system Galvanometer operation Styles ribbon Function of heat sensitive paper Electromotive force	
COMMUNICATIONS	
<u>PERFORMANCE MODES</u> Reading	<u>EXAMPLES</u> Operation manual
	<u>SKILLS/CONCEPTS</u> Comprehension

(TASK STATEMENT) OPERATE ELECTROCARDIOGRAPH

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY - HAZARD
<p>Patient Correctly set-up electrocardiograph</p>	<p>Set lead selector switch Record standardization Adjust baseline Record 6 limb leads Code each lead Record 6 chest leads</p>	<p><u>Safety</u> No X-ray or diathermy equipment must be operating near patient Check speed of paper occasionally</p> <p><u>Hazard</u> Electrical interference</p>
<p><u>DECISIONS</u> Identify electrical interference Identify tremor Determine proper hook-up Identify improper reading</p>	<p><u>CUES</u> No electrical inference should be seen (sawtooth notching) Somatic tremor (irregular peaks) caused by patient tenseness Wondering baseline caused by patient movement or drag on lead wires Jittery baseline caused by poor electrical contact Poorly defined baseline caused by stylus heat, contact with paper, insufficient pressure Standardization must be adjusted if voltage deflections exceed width of paper</p>	<p><u>ERRORS</u> Improper reading</p>

ASK STATEMENT) OPERATE ELECTROCARDIOGRAPH

SCIENCE	MATH -- NUMBER SYSTEMS																								
<p>Function and operation of heart Conduction system of heart - cardiac waves Anatomy of chest and heart Characteristic electrical pattern Positions of chest leads Marking system code</p>	<p>Coding system to mark leads</p> <table border="0"> <tr><td>I</td><td>.</td></tr> <tr><td>II</td><td>..</td></tr> <tr><td>III</td><td>...</td></tr> <tr><td>AVR</td><td>..-</td></tr> <tr><td>AVL</td><td>..-</td></tr> <tr><td>AVF</td><td>.....</td></tr> <tr><td>V₁</td><td>..</td></tr> <tr><td>V₂</td><td>..</td></tr> <tr><td>V₃</td><td>...</td></tr> <tr><td>V₄</td><td>...</td></tr> <tr><td>V₅</td><td>...</td></tr> <tr><td>V₆</td><td>...</td></tr> </table>	I	.	II	..	III	...	AVR	..-	AVL	..-	AVF	V ₁	..	V ₂	..	V ₃	...	V ₄	...	V ₅	...	V ₆	...
I	.																								
II	..																								
III	...																								
AVR	..-																								
AVL	..-																								
AVF																								
V ₁	..																								
V ₂	..																								
V ₃	...																								
V ₄	...																								
V ₅	...																								
V ₆	...																								
COMMUNICATIONS																									
<p><u>PERFORMANCE MODES</u></p> <p>Reading</p>	<p><u>EXAMPLES</u></p> <p>Machine operating manual, graph</p> <p><u>SKILLS/CONCEPTS</u></p> <p>Comprehension, description of mechanism, terminology</p>																								

(TASK STATEMENT) PERFORM FOLLOW-UP CARE OF PATIENT

244

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
Patient Towel Water	Remove cable cords Remove electrodes Wipe off jelly	Safety Clean off all jelly Hazard Jelly, although harmless and non- corrosive, is uncomfortable when left to dry on skin
<u>DECISIONS</u> Determine if patient needs extra attention	<u>CUES</u> Patient attitude Mental state	<u>ERRORS</u> Depression Fear

ASK STATEMENT) PERFORM FOLLOW-UP CARE OF PATIENT

SCIENCE	MATH - NUMBER SYSTEMS
COMMUNICATIONS	
<u>PERFORMANCE MODES</u> Speaking	<u>EXAMPLES</u> Reassure patient
<u>SKILLS/CONCEPTS</u> Appropriate diction, clarity of expression, persuasion, gestures, poise	

(TASK STATEMENT) PERFORM FOLLOW-UP CARE OF MACHINE

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Electrocardiograph Water Towels Pen Cleansing powder</p>	<p>Turn off power Remove tracing from machine Label tracing Clean electrodes with hot water and dry well Unplug power cord Cover machine</p>	<p>Safety Remove jelly from electrodes completely Hazard Tarnish and corrosion</p>
<p><u>DECISIONS</u> Select cleaning procedure</p>	<p><u>CUES</u> Use cleansing powder on electrodes Never steel wool</p>	<p><u>ERRORS</u> Improper operation</p>

244

(TASK STATEMENT) MOUNT ELECTROCARDIOGRAM

TOOLS, EQUIPMENT, MATERIALS, OBJECTS ACTED UPON	PERFORMANCE KNOWLEDGE	SAFETY -- HAZARD
<p>Properly labeled tracing Mounting folder</p>	<p>Select sections to be mounted Mount in appropriate spaces Label mounting folder with all essential information Include standardization in each lead</p>	
<p><u>DECISIONS</u></p> <p>Include abnormal beats, but not artifacts</p>	<p><u>CUES</u></p>	<p><u>ERRORS</u></p>

MOINT ELECTROCARDIOGRAM

ASK STATEMENT)

SCIENCE		MATH — NUMBER SYSTEMS	
Marking system code Characteristic electrical pattern of beats Characteristics of artifacts			
COMMUNICATIONS			
<u>PERFORMANCE MODES</u>		<u>EXAMPLES</u>	<u>SKILLS/CONCEPTS</u>
Writing		Label mounting folder name hospital number date room number time taken initials of technician doctor 245	Penmanship, spelling, accuracy, informational reports 14